

Datasheet: MCA6106

Description:	HUMAN ANTI IL6R (TOCILIZUMAB BIOSIMILAR)
Specificity:	IL6R
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	rhPM-1
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Human IgG1 (kappa) antibody - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Proclin 300
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Mouse AUK 12-20 antibody CDRs grafted onto human framework
External Database Links	UniProt: A0N0L5 Related reagents

Specificity

Human Anti-IL6R Antibody, clone rhPM-1 is a research grade biosimilar of the monoclonal antibody drug tocilizumab. It is a recombinant human IgG1 kappa antibody with variable regions from the therapeutic antibody tocilizumab.

Human Anti-IL6R Antibody is for research use only. It is suitable as a reference standard in a pharmacokinetic (PK) bridging ELISA with anti-tocilizumab antibodies, for example [HCA252](#) and [HCA253P](#). It can also be used as capture and detection reagent in an anti-drug antibody (ADA) bridging ELISA, with an anti-tocilizumab antibody in IgG1 format as reference standard, for example [HCA253](#).

Tocilizumab (Actemra/RoActemra) is a humanized IgG1/kappa monoclonal antibody approved for the treatment of moderate to severe rheumatoid arthritis and for the treatment of systemic juvenile idiopathic arthritis. The therapeutic antibody is directed against IL-6R and acts by blocking the binding of IL-6R (both the soluble form, sIL-6R, and the membrane bound form, mIL-6R) to IL-6, thus impeding the pro-inflammatory properties of this cytokine.

This product is NOT FOR THERAPEUTIC USE.

[View a summary of supporting anti-tocilizumab antibody products](#)

Purity	>98% by SDS PAGE
---------------	------------------

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10199 available at: 10199: https://www.bio-rad-antibodies.com/uploads/MSDS/10199.pdf
--------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Licensed Use	For in vitro research purposes only. Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
---------------------	--------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M340246:190104'

Printed on 11 Oct 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)