

## Datasheet: MCA6115APC

<b>Description:</b>	MOUSE ANTI HUMAN CD11a:APC
<b>Specificity:</b>	CD11a
<b>Other names:</b>	INTEGRIN ALPHA L CHAIN, LFA-1
<b>Format:</b>	APC
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	MEM-25
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 TESTS/1ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	Purified IgG conjugated to Allophycocyanin (APC) - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	APC	650	661
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative</b>	<0.1% Sodium Azide (NaN <sub>3</sub> )		
<b>Stabilisers</b>	0.2% Bovine Serum Albumin		
<b>Immunogen</b>	Leucocytes from a patient suffering from a LGL-type leukemia		
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P20701</a> <a href="#">Related reagents</a>		
	<b>Entrez Gene:</b> <a href="#">3683</a> ITGAL <a href="#">Related reagents</a>		

**Synonyms** CD11A

**Specificity** **Mouse anti Human CD11a antibody, clone MEM-25** recognizes the [integrin alphaL subunit](#), a 170-180 kDa type I transmembrane glycoprotein that associates with CD18 to form the leucocyte function-associated antigen-1, an alphaLbeta2 integrin expressed at the surface of all human leukocyte.

LFA-1 is a receptor for ICAMs, JAM-1 and F11R. Its interaction with endothelial ligands, such as ICAM-1 and JAM-1 is responsible for leukocyte arrest and extravasation ([Osterman et al.2002](#)). In addition to its role in leukocyte adhesion and transmigration, LFA-1 is involved in T cell differentiation ([Verma et al. 2016](#), [Verman and Kelleher 2017](#)) and is required for the effective lysis of infected cells by cytotoxic T cells ([Anikeeva et al. 2005](#), [Pontes Ferreira et al. 2017](#)). LFA-1 also contributes to NK cell cytotoxicity ([Domingo et al. 2004](#)) and participates to the clearance of apoptotic neutrophils by macrophages ([Kristóf et al. 2013](#)). The interaction of LFA-1 with ICAM-1 is the target of Lifitegrast, a drug approved for the treatment of Dry Eye Disease ([Abidi et al. 2016](#)).

Mouse anti Human CD11 antibody, clone MEM-25 inhibits the cytotoxic activity of NK cells ([Bazil et al. 1990](#)).

**Purity** >95% by SDS PAGE

**Flow Cytometry** Use 10ul of the undiluted product to label 1x10<sup>6</sup> cells in 100ul

**Storage** Store at +4°C. DO NOT FREEZE.  
This product should be stored undiluted. This product is photosensitive and should be protected from light.

**Guarantee** 12 months from date of despatch

**Health And Safety Information** Material Safety Datasheet documentation #10041 available at:  
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

**Regulatory** For research purposes only

## Related Products

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:APC \(MCA928APC\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M353558:190509'

Printed on 01 May 2020