

## Datasheet: MCA6121PE

<b>Description:</b>	MOUSE ANTI HUMAN CD150:RPE
<b>Specificity:</b>	CD150
<b>Other names:</b>	SLAM
<b>Format:</b>	RPE
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	SLAM.4
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 TESTS/1ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Target Species

Human

#### Product Form

Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid

#### Max Ex/Em

Fluorophore	Excitation Max (nm)	Emission Max (nm)
RPE 488nm laser	496	578
RPE 561nm laser	546	578

#### Preparation

Purified IgG prepared by ion exchange chromatography from tissue culture supernatant

#### Buffer Solution

Phosphate buffered saline

#### Preservative Stabilisers

<0.1% Sodium Azide (NaN<sub>3</sub>)  
0.2% Bovine Serum Albumin

#### Immunogen

Human CD150-transfected 300.19 cells

**External Database  
Links**

**UniProt:**

[Q13291](#)   [Related reagents](#)

**Entrez Gene:**

[6504](#)   SLAMF1   [Related reagents](#)

---

**Synonyms**

SLAM

---

**Specificity**

**Mouse anti Human CD150 antibody, clone SLAM.4** recognizes CD150 also known as signaling lymphocyte activation molecule (SLAM), is a member of the immunoglobulin superfamily. CD150 is expressed in dendritic cells, Natural Killer (NK) cells, B and T cells ([Bhattacharjee and Yadava 2018](#)).

CD150 is a membrane phosphoglycoprotein which contains 2 extracellular domains corresponding to a variable and constant region, a general feature of the Ig superfamily ([Cocks et al. 1995](#)). CD150 also has a cytoplasmic domain which contains 3 unique tyrosine motifs with a conserved TxYxxV sequence. This tyrosine motif is the binding domain for SLAM-associated protein (SAP) which form a complex together ([Sayos et al 1998](#), [Poy et al. 1999](#)). This complex plays a role in the regulation of Natural killer (NK) and T cell development. Mutations to SAP in humans results in X-linked lymphoproliferative disease which is a fatal immunological disorder ([Coffey et al. 1998](#)). Activation of CD150 co-activates T- or B-lymphocytes responses and plays a role in bidirectional T-B cell stimulation ([Carballido et al. 1997](#)). It can also regulate NK and T cell development, B cell regulation and antibody production/isotype switching ([Veillette 2006](#)). CD150 is also the receptor that the measles virus recognizes ([Tatsuo et al. 2000](#)).

The SLAM.4 clone has been used in flow cytometry to analyze the cell distribution of CD150 on leucocytes and lymphocytes as well as to distinguish between subsets of CD4+ and CD8+ T lymphocytes ([Romero et al. 2001](#)).

---

**Purity**

>95% by SDS PAGE

---

**Flow Cytometry**

Use 10ul of the undiluted reagent to label  $1 \times 10^6$  cells in 100ul

---

**Storage**

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light.

---

**Guarantee**

Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10041 available at:  
<https://www.bio-rad-antibodies.com/SDS/MCA6121PE>  
10041

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA928PE\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M401918:220718'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)