

Datasheet: MCA6156GA

BATCH NUMBER 157929

Description:	MOUSE ANTI HUMAN C5/C5a
Specificity:	C5/C5a
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	557
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/50
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
Immunogen	Human C5

**External Database
Links**

UniProt:

[P01031](#) [Related reagents](#)

Entrez Gene:

[727](#) C5 [Related reagents](#)

Synonyms

CPAMD4

Specificity

Mouse anti Human C5/C5a antibody, clone 557, recognizes an epitope on C5 and C5a. C5 is a component of all three complement activation pathways (classical, lectin and alternative), which lead to its cleavage by C5 convertase. Cleavage of C5 results in the formation of C5b, which initiates pathogenic lysis and C5a, a potent inflammatory mediator. C5a is a ligand for the G protein coupled receptors, C5aR1 and C5aR2 (also known as C5L2) through which its inflammatory functions are controlled ([Merle et al. 2015](#)). These receptors are expressed ubiquitously but are found at its highest levels on the surface of immune cells, including macrophages, neutrophils and T cells. The inflammatory responses mediated by C5a are chemotaxis, leukocyte degranulation, vascular permeability and cytokine plus chemokine production. Beyond innate immunity, C5a plays a role in organ development, tissue regeneration and hematopoiesis ([Zhang et al. 2017](#)).

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA6156GA>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC

Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M387833:210726'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)