

Datasheet: MCA6177B

BATCH NUMBER 100005275

Description:	MOUSE ANTI HUMAN CCL13:Biotin
Specificity:	CCL13
Other names:	MCP-4
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	B16-3D8
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide (NaN ₃)
Stabilisers	1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
Immunogen	<i>E. coli</i> derived recombinant human full length CCL13 (Phe17-Thr98)

External Database**Links****UniProt:**[Q99616](#)[Related reagents](#)**Entrez Gene:**[6357](#)

CCL13

[Related reagents](#)**Synonyms**

MCP4, NCC1, SCYA13

Fusion Partners

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

Specificity

Mouse anti Human CCL13 antibody, clone B16-3D8, recognizes CCL13, a member of the CC family of cytokines. CCL13, also known as MCP-4, is a 98 amino acid (aa) precursor protein which contains a 23 aa hydrophobic signal peptide that is cleaved to the mature form of the protein which is 75aa protein (8kDa). CCL13 is a chemoattractant for various cells including monocytes, T cells, immature dendritic cells and eosinophils ([Entiquez & Garcia-Zepeda 2013](#)). It has also been associated with various inflammatory diseases including asthma, arthritis, atherosclerosis and multiple sclerosis ([Procaccini *et al.* 2015](#)). CCL13 can bind to several chemokine receptors (CCR1, CCR2 and CCR3) to elicit chemoattractant responses. It is also known to induce degranulation of eosinophils and the secretion of pro-inflammatory cytokines.

This biotinylated Mouse anti Human CCL13 antibody, clone B16-3D8 (MCA6177B) can be used as a detection antibody in a sandwich ELISA with the purified Mouse anti Human CCL13 antibody, clone B06-1G9 ([MCA6176GA](#)) as the capture antibody.

Storage

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10041 available at:

<https://www.bio-rad-antibodies.com/SDS/MCA6177B>

10041

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE ANTI HUMAN CCL13 \(MCA6176GA\)](#)**North & South America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M378776:210304'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)