

Datasheet: MCA6184GA

Description:	MOUSE ANTI HUMAN CCL27
Specificity:	CCL27
Other names:	CTACK
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	A01-2H4
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen *E. coli* derived recombinant human full length CCL27 (Phe25-Gly112)

External Database

Links

UniProt:

[Q9Y4X3](#)

[Related reagents](#)

Entrez Gene:

[10850](#)

CCL27

[Related reagents](#)

Synonyms

ILC, SCYA27

Fusion Partners

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

Specificity

Mouse anti Human CCL27 antibody, clone A01-2H4, recognizes C-C motif chemokine 27, also known as CTACK, ESkinE, ILC, and SCYA27. CCL27 is a C-C motif chemokine. CCR10 is the only known receptor for CCL27, and is expressed in normal skin ([Mao et al. 2018](#)). CCL27 is able to attract a subset of memory / effector T cells from peripheral blood to specific locations on the skin ([Wang et al. 2018](#)). CCL27 plays a role in T cell mediated skin inflammation, and high levels have been found in patients with psoriasis and atopic dermatitis. Upregulation of CCL27 has been identified in patients with myeloma, and is correlated with poor overall survival ([Thangavadivel et al. 2016](#)). CCL27 is secreted by primary tumor cells with high metastatic potential and is able to attract lymphatic endothelial cells, a crucial process during lymphangiogenesis ([Karnezis et al. 2019](#)).

This purified Mouse anti Human CCL27 antibody, clone A01-2H4 (MCA6184GA) can be used as a capture antibody in a sandwich ELISA with the biotinylated Mouse anti Human CCL27 antibody, clone B05-10H9 ([MCA6185B](#)) as the detection antibody.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)

[HRP](#)

Rabbit Anti Mouse IgG (STAR12...)

[RPE](#)

Rabbit Anti Mouse IgG (STAR8...)

[DyLight®800](#)

Goat Anti Mouse IgG (STAR76...) [RPE](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®680](#),
[DyLight®800](#), [FITC](#), [HRP](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Recommended Useful Reagents

[MOUSE ANTI HUMAN CCL27:Biotin \(MCA6185B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M384205:210513'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)