

Datasheet: MCA6190GA

BATCH NUMBER 100005290

Description:	MOUSE ANTI HUMAN FIBRINOGEN
Specificity:	FIBRINOGEN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AB05-1F11
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Fibrinogen purified from human serum

External Database**Links****UniProt:**[P02679](#) [Related reagents](#)[P02671](#) [Related reagents](#)[P02675](#) [Related reagents](#)**Entrez Gene:**[2266](#) FGG [Related reagents](#)[2243](#) FGA [Related reagents](#)[2244](#) FGB [Related reagents](#)

Fusion Partners

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

Specificity

Mouse anti Human fibrinogen antibody, clone AB05-1F11 recognizes fibrinogen, a 340 kDa glycoprotein secreted by the liver which circulates in plasma as a soluble homodimer. The identical monomer subunits are linked by 17 disulfide bond and each subunit is composed of 3 polypeptide chains referred to as alpha, beta and gamma ([Tiscia & Margaglione 2018](#)).

In response to bleeding, the coagulation cascade is activated which results in the conversion of fibrinogen to fibrin. The resulting fibrin spontaneously forms polymers creating an insoluble gel which combines with platelets to form a spongy mass. This mass is the blood clot which subsequently hardens to prevent further bleeding from the blood vessel ([Tiscia & Margaglione 2018](#)). The impairment of this conversion from fibrinogen to fibrin is associated with a range of pathologies including coagulopathies, ischemic stroke and obstetrical complications. It is also implicated in various neurological diseases such as Alzheimer's, multiple sclerosis and traumatic central nervous system injury ([Petersen et al. 2018](#)).

The purified Mouse anti Human fibrinogen antibody, clone AB05-1F11 (MCA6190GA) can be used as a capture antibody in a sandwich ELISA with the biotinylated Mouse anti Human fibrinogen antibody, clone AB10-3C12 (MCA6191B) as the detection antibody.

Storage

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA6190GA>
10040

RegulatoryFor research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

Recommended Useful Reagents

[MOUSE ANTI HUMAN FIBRINOGEN:Biotin \(MCA6191B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M378783:210304'

Printed on 18 Jan 2024