

Datasheet: MCA6203B

BATCH NUMBER 100005093

Description:	MOUSE ANTI HUMAN CCL20:Biotin
Specificity:	CCL20
Other names:	MIP-3a
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	O03-3H12
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide (NaN ₃)
Stabilisers	1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
Immunogen	<i>E. coli</i> derived recombinant human full length CCL20 (Ala27-Met96)

External Database**Links****UniProt:**

[P78556](#) [Related reagents](#)

Entrez Gene:

[6364](#) CCL20 [Related reagents](#)

Synonyms

LARC, MIP3A, SCYA20

Fusion Partners

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

Specificity

Mouse anti Human CCL20 antibody, clone O03-3H12, recognizes C-C motif chemokine 20, also known as LARC, MIP3A, and SCYA20. The chemokine CCL20 is expressed mostly in the lymphatic tissue, lung, and liver, and is produced by cells involved with the inflammatory and immune responses. CCL20 plays an important role in recruiting inflammatory cells to the site of infection or injury ([Soto et al. 2017](#)). However, CCL20 is involved with pathogenesis of chronic inflammatory disorders, specifically interstitial lung fibrosis and chronic obstructive pulmonary disease (COPD), and can contribute to tumor metastasis ([Zhang et al. 2017](#)). Inhibition of CCL20, along with its ligand CCR6, might have valuable therapeutic potential for treatment of autoimmune and inflammatory diseases, and is therefore a target of interest for drug discovery ([Ranasinghe and Eri 2018](#)).

The biotinylated Mouse anti Human CCL20 antibody, clone O03-3H12 (MCA6203B) can be used as a detection antibody in a sandwich ELISA with the purified Mouse anti Human CCL20 antibody, clone J03-4D4 ([MCA6202GA](#)) as the capture antibody.

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10041 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA6203B>
10041

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE ANTI HUMAN CCL20 \(MCA6202GA\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M378788:210304'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)