

Datasheet: MCA6205B

**BATCH NUMBER 100006918**

<b>Description:</b>	MOUSE ANTI HUMAN INTERLEUKIN-11:Biotin
<b>Specificity:</b>	IL-11
<b>Other names:</b>	AGIF
<b>Format:</b>	Biotin
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	C16-4B8
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG conjugated to Biotin - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Stabilisers</b>	1% Bovine Serum Albumin
<b>Approx. Protein Concentrations</b>	IgG concentration 0.1 mg/ml
<b>Immunogen</b>	<i>E. coli</i> derived recombinant human full length interleukin-11 (Pro22-Leu199)

---

**External Database****Links****UniProt:**[P20809](#)[Related reagents](#)**Entrez Gene:**[3589](#)

IL11

[Related reagents](#)

---

**Fusion Partners**

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

---

**Specificity**

**Mouse anti Human interleukin-11 antibody, clone C16-4B8**, recognizes interleukin-11 (IL-11), also known as AGIF. IL-11 belongs to the IL-6 family of cytokines, and plays a protective role in acute inflammatory settings ([Traber et al. 2019](#)). IL-11 has both pro- and anti-inflammatory functions through activation of signaling cascades ([Lokau et al. 2017](#)). IL-11 is largely absent in the body fluids of healthy people, but is upregulated in a series of diseases. IL-11 can be expressed in fibroblasts, and its inhibition has been suggested as a treatment for fibrotic disease ([Schafer et al. 2017](#)), including fibrotic cardiovascular diseases ([Ye et al. 2019](#)). A role of IL-11 has also been identified in development of inflammation-associated cancers, and its inhibition may be a therapeutic goal ([Putoczki et al. 2013](#)).

The biotinylated Mouse anti Human interleukin-11 antibody, clone C16-4B8 (MCA6205B) can be used as a detection antibody in a sandwich ELISA with the purified Mouse anti Human interleukin-11 antibody, clone D08-4C12 ([MCA6204GA](#)) as the capture antibody.

---

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA6205B>  
10041

---

**Regulatory**

For research purposes only

---

**Related Products****Recommended Useful Reagents**[MOUSE ANTI HUMAN INTERLEUKIN-11 \(MCA6204GA\)](#)[RECOMBINANT HUMAN INTERLEUKIN-11 \(PHP303\)](#)**ELISA Matched Pair - Capture Antibody**[MOUSE ANTI HUMAN INTERLEUKIN-11 \(MCA6204GA\)](#)

## ELISA Antigen - Positive Control

### RECOMBINANT HUMAN INTERLEUKIN-11 (PHP303)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M384568:210513'

**Printed on 16 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)