

Datasheet: MCA6222GA

Description:	MOUSE ANTI HUMAN CXCL16
Specificity:	CXCL16
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	B06-3H6
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	<i>E. coli</i> derived recombinant human full length CXCL16 (Asn49-Pro137)

**External Database
Links**

UniProt:

[Q9H2A7](#) [Related reagents](#)

Entrez Gene:

[58191](#) CXCL16 [Related reagents](#)

Synonyms

SCYB16, SRPSOX

Fusion Partners

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

Specificity

Mouse anti Human CXCL16 antibody, clone B06-3H6, recognizes C-X-C motif Chemokine 16 (CXCL16). Unlike most chemokines, CXCL16 is a membrane bound ligand anchored through transmembrane domains. Although a soluble form of CXCL16 is generated through proteolytic cleavage. Transmembrane CXCL16 is mainly expressed on macrophages, dendritic cells, B cells and various tumor cells. In this form, it acts as an adhesion molecule though the binding of its receptor CXCR6 expressed on other cells including T cells. Transmembrane CXCL16 also binds gram positive and gram negative bacteria, playing a role in innate and adaptive immune responses ([Shimaoka et al. 2003](#)).

Soluble CXCL16 has chemoattractant activity, and elicits this response through the receptor CXCR6. CXCL16 may also play a role in tumor progression and serve as a biomarker for the evaluation of some tumor type prognoses ([Richardsen et al. 2015](#), [Shibata et al. 2020](#), [Gooden et al. 2014](#)).

The purified Mouse anti Human CXCL16 antibody, clone B06-3H6 (**MCA6222GA**) can be used as a capture antibody in a sandwich ELISA with the biotinylated Mouse anti Human CXCL16 antibody, clone C19-5F9 ([MCA6223B](#)) as the detection antibody.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA6222GA>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)

[HRP](#)

Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR8...)	DyLight@800
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight@488 , DyLight@550 , DyLight@650 , DyLight@680 , DyLight@800 , FITC , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Recommended Useful Reagents

[MOUSE ANTI HUMAN CXCL16:Biotin \(MCA6223B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M384103:210513'

Printed on 12 Aug 2023