

Datasheet: MCA6225B

BATCH NUMBER 100005236

Description:	MOUSE ANTI HUMAN CXCL7:Biotin
Specificity:	CXCL7
Other names:	Platelet basic protein
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	C06-A43
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
Immunogen	<i>E. coli</i> derived recombinant human full length CXCL7 (Ala59-Asp128)

External Database Links	<p>UniProt: P02775 Related reagents</p> <p>Entrez Gene: 5473 PPBP Related reagents</p> <hr/> <p>Synonyms CTAP3, CXCL7, SCYB7, TGB1, THBGB1</p> <hr/> <p>Fusion Partners Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0</p> <hr/> <p>Specificity Mouse anti Human CXCL7 antibody, clone C11-5G10, recognizes C-X-C motif chemokine 7 (CXCL7), also known as Platelet Basic Protein, a platelet derived growth factor belonging to the C-X-C subfamily of chemokines. Expression of CXCL7 is largely restricted to platelets where it is produced as an inactive precursor. Upon tissue injury or infection, CXCL7 is one of the first chemokines to appear at the focal sites to mediate neutrophil recruitment after it has been modified by a series of cleavage and activation steps. CXCL7 mediates neutrophil recruitment through binding of its receptors CXCR1 or CXCR2.</p> <p>CXCL7 has been noted as playing important roles in tumorigenesis and has been associated with the proliferation and metastasis of a range of tumor types (Brown et al. 2017). As such, there have been a range of studies performed testing the use of serum CXCL7 as a diagnostic biomarker in colorectal cancer, lung cancer and renal cell carcinoma samples (Li et al. 2019, Kinouchi et al. 2017, Du et al. 2018). These studies have suggested that serum CXCL7 may be a sensitive and specific biomarker for the diagnosis of these cancer types.</p> <p>The biotinylated Mouse anti Human CXCL7 antibody, clone C06-4A3 (MCA6225B) can be used as a detection antibody in a sandwich ELISA with the purified Mouse anti Human CXCL7 antibody, clone C11-5G10 (MCA6224GA) as the capture antibody.</p> <hr/> <p>Storage Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p> <hr/> <p>Guarantee 12 months from date of despatch</p> <hr/> <p>Health And Safety Information Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf</p> <hr/> <p>Regulatory For research purposes only</p>
--------------------------------	--

Related Products

[ELISA Matched Pair - Capture Antibody](#)

[MOUSE ANTI HUMAN CXCL7 \(MCA6224GA\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M386383:210521'

Printed on 08 Feb 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)