

Datasheet: MCA6244GA

Description:	MOUSE ANTI HUMAN INTERLEUKIN-18BP
Specificity:	IL-18BP
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	B01-3E12
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	NS0 derived recombinant human full length interleukin-18BP (Thr31-Gly194)

**External Database
Links**

UniProt:

[O95998](#) [Related reagents](#)

Entrez Gene:

[10068](#) IL18BP [Related reagents](#)

Fusion Partners

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

Specificity

Mouse anti Human interleukin-18BP antibody, clone B01-3E12, recognizes interleukin-18-binding protein (IL-18BP), also known as takedin- α . IL-18BP is a secreted immunoregulatory binding protein which binds to and inhibits the cytokine activity of mature IL-18 with high affinity. IL-18 plays a role in host defense, inflammation, tissue repair, and regeneration ([Mühl & Bachmann 2019](#)).

IL-18BP is constitutively expressed in the spleen, PBMCs, monocytes and umbilical vein endothelial cells ([Mühl & Bachmann 2019](#), [Corbaz *et al.* 2002](#)). It is detectable in the serum of healthy humans at a 20 fold molar excess compared to IL-18, with subjects showing levels of 2,000 – 4,000 pg/ml of IL-18BP ([Kaplanski *et al.* 2018](#)). IL-18BP is upregulated in response to inflammation as part of the inflammation resolution process. Additionally, it can be upregulated during various cancers and IL-18BP levels are also affected by viral infections such as dengue and hepatitis C ([Mühl & Bachmann 2019](#)).

The purified Mouse anti Human interleukin-18BP antibody, clone B01-3E12 (MCA6244GA) can be used as a capture antibody in a sandwich ELISA with the biotinylated Mouse anti Human interleukin-18BP antibody, clone A07-5D2 ([MCA6245B](#)) as the detection antibody.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...) [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)

Goat Anti Mouse IgG (STAR76...) [RPE](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®680](#),
[DyLight®800](#), [FITC](#), [HRP](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Recommended Useful Reagents

[MOUSE ANTI HUMAN INTERLEUKIN-18BP:Biotin \(MCA6245B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M384125:210513'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)