

Datasheet: MCA6268GA

BATCH NUMBER 100005443

Description:	MOUSE ANTI HUMAN SAP
Specificity:	SAP
Other names:	Serum amyloid P-component
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AB18-4E2
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen	Native SAP purified from human serum
External Database Links	<p>UniProt: P02743 Related reagents</p> <p>Entrez Gene: 325 APCS Related reagents</p>
Synonyms	PTX2
Fusion Partners	Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0
Specificity	<p>Mouse anti Human SAP antibody, clone AB18-4E2, recognizes Serum Amyloid P-component (SAP) also known as 9.5S alpha-1-glycoprotein, a member of the pentraxin family of proteins. It is synthesized by hepatocytes and secreted into the blood. In humans the levels of SAP range between 20-50 µg/ml (Pilling & Gomer 2018).</p> <p>SAP is able to bind to a variety of molecules including apoptotic debris, bacterial polysaccharides, amyloid deposits and bacterial toxins. Phagocytic cells such as monocytes and macrophages then bind the SAP and engulf the material bound to SAP as part of debris clearance (Pilling & Gomer 2018).</p> <p>The purified Mouse anti Human SAP antibody, clone AB18-4E2 (MCA62648GA) can be used as a capture antibody in a sandwich ELISA with the biotinylated Mouse anti Human SAP antibody, clone AB12-5A1 (MCA6269B) as the detection antibody.</p>
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA6268GA 10040</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP

Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Recommended Useful Reagents

[MOUSE ANTI HUMAN SAP:Biotin \(MCA6269B\)](#)

ELISA Matched Pair - Detection Antibody

[MOUSE ANTI HUMAN SAP:Biotin \(MCA6269B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M386398:210521'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)