

Datasheet: MCA6280GA

Description:	MOUSE ANTI HUMAN CD258
Specificity:	CD258
Other names:	TNFSF14
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	D04-8A4
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen	<i>E. coli</i> derived recombinant human CD258 (Leu61-Val240)
External Database Links	<p>UniProt: O43557 Related reagents</p> <p>Entrez Gene: 8740 TNFSF14 Related reagents</p>
Synonyms	HVEML, LIGHT
Fusion Partners	Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0
Specificity	<p>Mouse anti Human CD258 antibody, clone D04-8A4, recognizes CD258 also known as TNFSF14, Herpes virus entry mediator and LIGHT. CD258 is the ligand of the receptor TNFRSF14, which is the cell surface receptor that the Herpes virus uses to attach to the cell surface and to enter into the cell, it is also a ligand of lymphotoxin-β receptor (Ware & Sedý 2011). Signaling events initiated upon binding of CD258 to these receptors is cell dependent. However, the responses are often related to immune system stimulation including the production of various cytokines, adhesion molecules that control the migration and positioning of dendritic cells and lymphocytes (Skeate et al. 2020, Tang et al. 2016).</p> <p>CD258 is primarily expressed by activated T cells, activated natural killer cells and immature dendritic cells. It is initially synthesized as a type II transmembrane protein, which can be proteolytically cleaved to release a soluble bioactive form (Ware & Sedý 2011).</p> <p>The purified Mouse anti Human CD258 antibody, clone D04-8A4 (MCA62680GA) can be used as a capture antibody in a sandwich ELISA with the biotinylated Mouse anti Human CD258 antibody, clone B01-9B3 (MCA6281B) as the detection antibody.</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR76...)	RPE
Human Anti Mouse IgG2a (HCA037...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®680 , DyLight®800 , FITC , HRP
Goat Anti Mouse IgG (STAR70...)	FITC

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

Recommended Useful Reagents

[MOUSE ANTI HUMAN CD258:Biotin \(MCA6281B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M384559:210513'

Printed on 21 Mar 2022