

Datasheet: MCA6281B

Description:	MOUSE ANTI HUMAN CD258:Biotin
Specificity:	CD258
Other names:	TNFSF14
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	B01-9B3
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
Immunogen	<i>E. coli</i> derived recombinant human CD258 (Leu61-Val240)

**External Database
Links**

UniProt:

[O43557](#) [Related reagents](#)

Entrez Gene:

[8740](#) TNFSF14 [Related reagents](#)

Synonyms

HVEML, LIGHT

Fusion Partners

Cell fusion between immunized BALB/c mouse spleen cells and mouse myeloma SP2/0

Specificity

Mouse anti Human CD258 antibody, clone B01-9B3, recognizes CD258 also known as TNFSF14, Herpes virus entry mediator and LIGHT. CD258 is the ligand of the receptor TNFRSF14, which is the cell surface receptor that the Herpes virus uses to attach to the cell surface and to enter into the cell, it is also a ligand of lymphotoxin- β receptor ([Ware & Sedý 2011](#)). Signaling events initiated upon binding of CD258 to these receptors is cell dependent. However, the responses are often related to immune system stimulation including the production of various cytokines, adhesion molecules that control the migration and positioning of dendritic cells and lymphocytes ([Skeate et al. 2020](#), [Tang et al. 2016](#)).

CD258 is primarily expressed by activated T cells, activated natural killer cells and immature dendritic cells. It is initially synthesized as a type II transmembrane protein, which can be proteolytically cleaved to release a soluble bioactive form ([Ware & Sedý 2011](#)).

The biotinylated Mouse anti Human CD258 antibody, clone B01-9B3 (MCA62781B) can be used as a detection antibody in a sandwich ELISA with the purified Mouse anti Human CD258 antibody, clone D04-8A4 ([MCA6280GA](#)) as the capture antibody.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10041 available at:
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE ANTI HUMAN CD258 \(MCA6280GA\)](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M384128:210513'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)