

Datasheet: MCA6296

Description:	RABBIT ANTI HISTONE H3F3A (pThr11)
Specificity:	HISTONE H3F3A (pThr11)
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	RM164
Isotype:	IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1 µg/ml - 0.2 µg/ml
Western Blotting	▪			2 µg/ml - 0.5 µg/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Based on sequence similarity, is expected to react with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Purified IgG prepared by affinity chromatography of Protein A from animal origin-free culture supernatant

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide
1% Bovine Serum Albumin

50% Glycerol

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
---------------------------------------	-----------------------------

Immunogen	A phospho-peptide corresponding to H3F3A (pThr11)
------------------	---

External Database Links	UniProt: P84243 Related reagents
	Entrez Gene: 3021 H3F3B Related reagents

Synonyms	H3.3A, H3.3B, H3F3
-----------------	--------------------

Specificity	Rabbit anti Human histone H3F3A (pThr11) antibody recognizes histone H3.3, when phosphorylated at threonine 11. In eukaryotic cells, histones are proteins which associate with nuclear DNA and help condense it into chromatin, playing important roles in gene regulation. H3F3A is a common non-centromeric variant of histone H3, and appears to play a role in both gene activation and silencing (Shi et al. 2017). Post-translational modifications of histones alter the interactions between histones and DNA, and between histones and other proteins, with dramatic consequences for DNA transcription and other biological events (Sawika and Seiser 2012). This protein can undergo massive phosphorylation at threonine 11, among other sites, altering the function of histone H3 in a temporal and spatial way. This antibody detects H3F3A when phosphorylated at threonine 11.
--------------------	---

Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10048 available at: 10048: https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (STAR35...)	RPE
Sheep Anti Rabbit IgG (STAR36...)	DyLight@488 , DyLight@680 , DyLight@800

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M376276:210127'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)