

Datasheet: MCA6321

Description:	RABBIT ANTI CYCLIN B1
Specificity:	CYCLIN B1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	RM281
Isotype:	IgG
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			1/1000 - 1/2000
Western Blotting	▪			1/1000 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography of Protein A from animal origin-free culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Stabilisers	1% Bovine Serum Albumin 50% Glycerol
Immunogen	A peptide corresponding to cyclin B1
External Database Links	UniProt:

Entrez Gene:

[891](#) CCNB1 [Related reagents](#)

Synonyms	CCNB
-----------------	------

Specificity	<p>Rabbit anti Human cyclin B1 antibody recognizes G2/mitotic-specific cyclin-B1, also known as CCNB1. Cyclins are protein kinases required for cell cycle progression, and different cyclins execute their functions at different stages of the cell cycle by binding to cyclin-dependent kinases (CDKs). Cyclin B1 forms a complex with CDK1 and regulates entry into mitosis at the G2/M transition, synchronizing critical events such as nuclear envelope breakdown and centrosome separation (Wang et al.. 2014). Cyclin B1 and CDK1 are present in the mitochondria of mammalian cells, and their expression increases under conditions of DNA damage (Xie et al.. 2019). Overexpression of cyclin B1 in solid tumors has been associated with worse prognosis of lung cancer and esophageal cancer, but better prognosis of colorectal cancer (Ye et al.. 2016). Cyclin B1 interacts with BRCA1, which when mutated, carries a substantial risk of breast and ovarian cancer. Targeting Cyclin B1 has been suggested for BRCA1-associated breast cancer therapy (Choi et al.. 2018).</p>
--------------------	---

Storage	<p>Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use</p>
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	<p>Material Safety Datasheet documentation #10048 available at: https://www.bio-rad-antibodies.com/SDS/MCA6321 10048</p>
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M386190:210519'

