

Datasheet: MCA6337Z

**BATCH NUMBER 165243**

<b>Description:</b>	MOUSE ANTI HUMAN CCL11:Preservative Free
<b>Specificity:</b>	CCL11
<b>Other names:</b>	Eotaxin
<b>Format:</b>	Preservative Free
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	1.3_17H09-02C06
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	<b>Yes</b>	<b>No</b>	<b>Not Determined</b>	<b>Suggested Dilution</b>
Flow Cytometry	▪			
ELISA	▪			2.0 ug/ml - 1.0 ug/ml
Western Blotting	▪			0.25 ug/ml - 0.5 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - lyophilized
<b>Reconstitution</b>	Reconstitute with 0.5 ml distilled water
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present

<b>Carrier Free</b>	Yes
<b>Endotoxin Level</b>	< 1 EU/μg
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml after reconstitution
<b>Immunogen</b>	<i>E. coli</i> derived recombinant human CCL11
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P51671</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">6356</a>    CCL11    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	SCYA11
<b>Specificity</b>	<p><b>Mouse anti Human CCL11 antibody, clone 1.3_17H09-02C06</b>, recognizes C-C motif chemokine 11 (CCL11), also known as eotaxin, eosinophil chemotactic protein and small inducible cytokine A11. CCL11 is produced by a range of cell types including epithelial cells, fibroblasts, smooth muscle cells, astrocytes, chondrocytes and tissue resident macrophages. Additionally, in the central nervous system, it is expressed by choroid plexus epithelial cells, pericytes, and microglia under inflammatory conditions (<a href="#">Teixeira et al. 2018</a>).</p> <p>CCL11 acts as a chemoattractant for eosinophil recruitment to inflammatory sites, including during asthma attacks (<a href="#">Wen &amp; Rothenberg 2016</a>). Upon binding of CCL11 to its receptor CCR3, a series of signaling events are triggered which can result in calcium mobilization, CD11b upregulation, MAPK activation and actin polymerization. It can also result in a rapid change of shape associated with chemotaxis and granule release (<a href="#">Conroy &amp; Williams 2001</a>).</p>
<b>Storage</b>	<p>Prior to reconstitution store at -20°C.  After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10294 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA6337Z">https://www.bio-rad-antibodies.com/SDS/MCA6337Z</a> 10294
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR8...)	<a href="#">DyLight@800</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">Alk. Phos.</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight@488</a> , <a href="#">DyLight@680</a> , <a href="#">DyLight@800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M386539:210527'

Printed on 19 Jan 2024