

Datasheet: MCA6351

Description:	RABBIT ANTI SATB2
Specificity:	SATB2
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	RM365
Isotype:	IgG
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			1/200 - 1/1000
Western Blotting	▪			1/1000 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from animal origin-free culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Stabilisers	1% Bovine Serum Albumin 50% Glycerol
Immunogen	A peptide corresponding to the C-terminus of human DNA-binding protein SATB2
External Database Links	UniProt:

Entrez Gene:

[23314](#) SATB2 [Related reagents](#)

Synonyms KIAA1034

Specificity **Rabbit anti SATB2 antibody** recognizes DNA-binding protein SATB2, also known as special AT-rich sequence-binding protein 2.

SATB2 is a transcription factor which can bind to matrix-attachment regions (MARs) and activate transcription of multiple genes simultaneously. This allows SATB2 to act as a high-level regulator of multiple networks involved with development ([Zarate and Fish 2016](#)). SATB2 is highly expressed in the brain and plays a role in development of the central nervous system and neocortical organization. It also plays a role in differentiation of osteoblasts ([Berg et al. 2016](#)). Mutations in the gene encoding SATB2 can lead to SATB2-associated syndrome, characterized by developmental delay, intellectual disability, speech and behavioural problems, and craniofacial abnormalities ([Whitton et al. 2018](#)).

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10048 available at:
10048: <https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M386980:210604'

Printed on 07 Feb 2023
