

Datasheet: MCA6374

Description:	MOUSE ANTI HUMAN RNA POL II CTD (pSer2/pSer5)
Specificity:	RNA POL II CTD (pSer2/pSer5)
Other names:	RPB1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	1A12G10
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			0.5 µl/ml - 2 µl/ml
Immunofluorescence	▪			5 µg/ml
Chromatin Immunoprecipitation	▪			5 µg per ChIP - 10 µg per ChIP

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.035% Sodium Azide (NaN ₃) 30% Glycerol
Carrier Free	Yes

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	A synthetic peptide containing phosphorylated serine 2 and phosphorylated serine 5 of the CTD of human RNA polymerase II
External Database Links	<p>UniProt: P24928 Related reagents</p> <p>Entrez Gene: 5430 POLR2A Related reagents</p>
Synonyms	POLR2
Specificity	<p>Mouse anti Human RNA pol II CTD (pSer2/pSer5) antibody, clone 1A12G10 recognizes RNA polymerase II, also known as RNAP II and pol II, when phosphorylated at serine 2 or serine 5 of the carboxyl-terminal repeat domain.</p> <p>RNA pol II is a 550 kDa enzyme complex made up of 12 subunits that transcribes DNA into precursors of messenger RNA (mRNA), snRNA, and microRNA (Kornburg, 1999) (Sims et al. 2004). The carboxyl-terminal repeat domain (CTD) of RNA pol II is reversibly phosphorylated in during the transcription cycle (Phatnani and Greenleaf, 2006). Repeats phosphorylated on both serine 2 and serine 5 are required for optimal binding of the SRI domain of the methyl transferase SET2 to RNA pol II (Rebehmed et al. 2014).</p> <p>Wide species cross-reactivity is expected from Mouse anti Human RNA pol II CTD (pSer2/pSer5), clone 1A12G10 based on sequence.</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10049 available at: https://www.bio-rad-antibodies.com/SDS/MCA6374 10049
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M394790:220218'

Printed on 18 Jan 2024