

Datasheet: MCA686S

BATCH NUMBER 170505

| | |
|----------------------|--------------------------------|
| Description: | MOUSE ANTI HUMAN MBP (aa84-89) |
| Specificity: | MBP (aa84-89) |
| Other names: | MYELIN BASIC PROTEIN |
| Format: | S/N |
| Product Type: | Monoclonal Antibody |
| Clone: | 22 |
| Isotype: | IgG2b |
| Quantity: | 2 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | ▪ | | | |
| Immunohistology - Paraffin | ▪ | | | |
| ELISA | | | ▪ | |
| Immunoprecipitation | | | ▪ | |
| Western Blotting | | | ▪ | |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|--|
| Target Species | Human |
| Product Form | Tissue Culture Supernatant - liquid |
| Preservative Stabilisers | 0.05% Sodium Azide |
| Immunogen | Human myelin basic protein. |
| External Database Links | UniProt: P02686 Related reagents |

Entrez Gene:[4155](#) MBP [Related reagents](#)

| | |
|-----------------------------------|--|
| RRID | AB_325006 |
| Fusion Partners | Spleen cells from immunized mice were fused with cells of the mouse SP2/0 myeloma cell line. |
| Specificity | Mouse anti MBP antibody, clone 22 recognizes an epitope between amino acids 84 and 89 of human myelin basic protein, also known as MBP, Myelin A1 protein or Myelin membrane encephalitogenic protein. MBP is a 304 amino acid ~33 kDa abundant component of the myelin membrane. There are 6 potential isoforms of MBP generated by alternative splicing, 4 of which lack the N-terminal region containing the immunogen sequence. Mouse anti MBP antibody, clone 22 is expected to recognize isoforms 1 and 2 (UniProt: P02686). The N-terminal region of MBP is highly conserved and Mouse anti MBP antibody, clone 22 is expected to be broadly species cross reactive for isoforms maintaining the N-terminal region. Mouse anti MBP antibody, clone 22 has been reported as suitable for use in western blotting (Chignola et al. 2000). |
| Histology Positive Control Tissue | Brain |
| References | <ol style="list-style-type: none">1. Groome, N. <i>et al.</i> (1988) New monoclonal antibodies reactive with defined sequential epitopes in human myelin basic protein. J Neuroimmunol. 19 (4): 305-15.2. Martenson, R. E. (1984) Experimental allergic encephalomyelitis. A useful model of multiple sclerosis. Eds Alvord, Kies and Suckling pp 511-521. Alan Liss N.Y.3. Harrower, T.P. <i>et al.</i> (2006) Long-term survival and integration of porcine expanded neural precursor cell grafts in a rat model of Parkinson's disease. Exp Neurol. 197 (1): 56-69.4. Liao, M.C. <i>et al.</i> (2009) Degradation of amyloid beta protein by purified myelin basic protein. J Biol Chem. 284: 28917-25.5. Diemel, L.T. <i>et al.</i> (2003) Role for TGF-beta1, FGF-2 and PDGF-AA in a myelination of CNS aggregate cultures enriched with macrophages. J Neurosci Res. 74 (6): 858-67. |
| Storage | <p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p> |
| Guarantee | 12 months from date of despatch |
| Health And Safety Information | Material Safety Datasheet documentation #10055 available at: https://www.bio-rad-antibodies.com/SDS/MCA686S 10055 |
| Regulatory | For research purposes only |

Related Products

Recommended Secondary Antibodies

| | |
|---|---|
| Rabbit Anti Mouse IgG (STAR12...) | RPE |
| Goat Anti Mouse IgG IgA IgM (STAR87...) | HRP |
| Goat Anti Mouse IgG (STAR76...) | RPE |
| Rabbit Anti Mouse IgG (STAR13...) | HRP |
| Goat Anti Mouse IgG (STAR70...) | FITC |
| Goat Anti Mouse IgG (H/L) (STAR117...) | Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP |
| Rabbit Anti Mouse IgG (STAR9...) | FITC |
| Goat Anti Mouse IgG (STAR77...) | HRP |
| Goat Anti Mouse IgG (Fc) (STAR120...) | FITC , HRP |

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M368765:200529'

Printed on 18 Jan 2024