

Datasheet: MCA70

BATCH NUMBER 155730

Description:	MOUSE ANTI BOVINE MBP (aa129-138)
Specificity:	MBP (aa129-138)
Other names:	MYELIN BASIC PROTEIN
Format:	S/N
Product Type:	Monoclonal Antibody
Clone:	1
Isotype:	IgG2a
Quantity:	2 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			1/10
Immunohistology - Paraffin			▪	
ELISA	▪			1/200 - 1/1000
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Bovine

Species Cross Reactivity

Reacts with: Rat, Human

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Tissue Culture Supernatant - liquid

Preparation

Tissue Culture Supernatant containing 0.1M Tris/HCl pH7.4 and 5-10% foetal calf serum

Preservative Stabilisers	0.05% Sodium Azide
Immunogen	Bovine Myelin Basic Protein (MBP).
External Database Links	<p>UniProt: P02687 Related reagents</p> <p>Entrez Gene: 618684 MBP Related reagents</p>
RRID	AB_2140358
Fusion Partners	Spleen cells from immunised NIH mice were fused with cells of the NS0 mouse myeloma cell line.
Specificity	<p>Mouse anti Bovine MBP (aa129-138) antibody, clone 1 recognizes myelin basic protein (MBP). Clone 1 is reactive with an epitope in the 129-138 region of the human MBP molecule. The numbering of MBP residues is that as described by Martenson Martenson, R. E., 1984.</p> <p>Clone 1 (129-138) has been reported as being suitable for Western blotting.</p>
References	<ol style="list-style-type: none"> 1. Groome, N. <i>et al.</i> (1985) Preparation and properties of monoclonal antibodies to myelin basic protein and its peptides. Neurochem Int. 7 (2): 309-17. 2. Martenson, R. E. (1984) Experimental allergic encephalomyelitis. A useful model for multiple sclerosis. In Alvord, E. C. <i>et al.</i> (Eds). Wiley, New York, pp273-289. 3. Groome, N.P. <i>et al.</i> (1986) Region-specific immunoassays for human myelin basic protein. J Neuroimmunol. 12 (4): 253-64. 4. Marignier, R. <i>et al.</i> (2010) Oligodendrocytes are damaged by neuromyelitis optica immunoglobulin G via astrocyte injury. Brain. 133 (9): 2578-91. 5. Brunner, C. <i>et al.</i> (1989) Differential ultrastructural localization of myelin basic protein, myelin/oligodendroglial glycoprotein, and 2',3'-cyclic nucleotide 3'-phosphodiesterase in the CNS of adult rats. J Neurochem. 52: 296-304. 6. Jatana, M. <i>et al.</i> (2006) Combination of systemic hypothermia and N-acetylcysteine attenuates hypoxic-ischemic brain injury in neonatal rats. Pediatr Res. 59: 684-9. 7. Chujor, C.S. <i>et al.</i> (1991) Serum antibodies against peripheral nervous system antigens in leprosy. Int J Lepr Other Mycobact Dis. 59: 590-7.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10336 available at: https://www.bio-rad-antibodies.com/SDS/MCA70 10336
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M374832:201201'

Printed on 19 Oct 2023