

Datasheet: MCA752F

BATCH NUMBER 171657

Description:	MOUSE ANTI GUINEA PIG CD8:FITC
Specificity:	CD8
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	CT6
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Guinea Pig		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin		

Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
Immunogen	Guinea pig peritoneal T-cells.
RRID	AB_321399
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the X63.Ag8.653 mouse myeloma cell line.
Specificity	Mouse anti Guinea Pig CD8 antibody, clone CT6 reacts with guinea pig CD8 present on cytotoxic T-cells. CD8 comprises 2 subunits, alpha and beta and exists as either an alpha/alpha homodimer or an alpha/beta heterodimer. Sequence analysis suggests that guinea pig CD8 is more closely related to human than rat or mouse CD8 (Nagarajan <i>et al.</i> 2004).
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul
References	<ol style="list-style-type: none"> 1. Tan, B.T.G. <i>et al.</i> (1985) Production of monoclonal antibodies defining guinea pig T-cell surface markers and a strain 13 Ia-like antigen: the value of immunohistological screening. Hybridoma 4: 115-124. 2. Steerenberg, P.A. <i>et al.</i> (1991) Tumour rejection after transfer of line 10 immunity is mediated by two T-cell populations. Cancer Immunol. Immunother. 34: 103-110. 3. Steerenberg, P.A. <i>et al.</i> (1990) Tumour infiltrating leucocytes (tils) during progressive tumour growth and BCG - mediated tumour regression. Virchows Archiv Cell Pathol. 59: 185-194. 4. Baker, D. <i>et al.</i> (1987) Changes in lymphocyte subsets after treatment with cyclophosphamide and during the development of contact sensitivity in the guinea pig. Int J Immunopharmacol. 9 (2): 175-83. 5. Antoniou, A.V. <i>et al.</i> (1986) Immunocytochemical identification and quantitation of mononuclear cells in the meninges during the development of chronic relapsing experimental allergic encephalomyelitis (CREAE) in the guinea pig. Cell Immunol. 97 (2): 386-96. 6. Liversidge, J. <i>et al.</i> (1987) EAU in the guinea pig: inhibition of cell-mediated immunity and Ia antigen expression by cyclosporin A. Clin Exp Immunol. 69 (3): 591-600. 7. Liversidge, J. & Forrester, J.V. (1988) Experimental autoimmune uveitis (EAU): immunophenotypic analysis of inflammatory cells in chorio retinal lesions. Curr Eye Res. 7 (12): 1231-41. 8. Debout, C. <i>et al.</i> (1991) The Kurloff cell in estrogenized guinea pigs as a CT7+ 8BE6-CT6- MR-1- CT10- IgM- lymphocyte with natural killer activity. Nat Immun Cell Growth Regul. 10 (6): 327-35. 9. Shang S <i>et al.</i> (2011) Activities of TMC207, rifampin, and pyrazinamide against Mycobacterium tuberculosis infection in guinea pigs. Antimicrob Agents Chemother. 55 (1): 124-31. 10. Lacy HM <i>et al.</i> (2011) Essential role for neutrophils in pathogenesis and adaptive immunity in <i>Chlamydia caviae</i> ocular infections. Infect Immun. 79 (5): 1889-97. 11. Komori, T. <i>et al.</i> (2011) A Microbial Glycolipid Functions as a New Class of Target Antigen for Delayed-type Hypersensitivity. J Biol Chem. 286: 16800-6.

12. Hiromatsu, K. *et al.* (2002) Induction of CD1-restricted immune responses in guinea pigs by immunization with mycobacterial lipid antigens. [J Immunol. 169: 330-9.](#)
13. Kramp, J.C. *et al.* (2011) The *in vivo* immunomodulatory effect of recombinant tumour necrosis factor-alpha in guinea pigs vaccinated with *Mycobacterium bovis* bacille Calmette-Guérin. [Clin Exp Immunol. 165: 110-20.](#)
14. Mishra, N.C. *et al.* (2010) Sulfur mustard induces immune sensitization in hairless guinea pigs. [Int Immunopharmacol. 10: 193-9.](#)
15. Schleiss, M.R. *et al.* (2007) Preconceptual administration of an alphavirus replicon UL83 (pp65 homolog) vaccine induces humoral and cellular immunity and improves pregnancy outcome in the guinea pig model of congenital cytomegalovirus infection. [J Infect Dis. 195: 789-98.](#)
16. Turner, O.C. *et al.* (2003) Immunopathogenesis of pulmonary granulomas in the guinea pig after infection with *Mycobacterium tuberculosis*. [Infect Immun. 71: 864-71.](#)
17. Wang, Y. *et al.* (2010) Local host response to chlamydial urethral infection in male guinea pigs. [Infect Immun. 78: 1670-81.](#)
18. Yang, H. *et al.* (2011) Three protein cocktails mediate delayed-type hypersensitivity responses indistinguishable from that elicited by purified protein derivative in the guinea pig model of *Mycobacterium tuberculosis* infection. [Infect Immun. 79: 716-23.](#)
19. Baldwin SL *et al.* (2012) The importance of adjuvant formulation in the development of a tuberculosis vaccine. [J Immunol. 188 \(5\): 2189-97.](#)
20. Obregón-Henao, A. Shang, s. Shanley, C.A. Basaraba, R.J. Caraway, M.L. Duncan, C.G. Ordway, D.J. Orme, I.M. (2013) Cortisone-Forced Reactivation of Weakly Acid Fast Positive Mycobacterium Tuberculosis in Guinea Pigs Previously Treated With Chemotherapy [Mycobacterial Diseases. 2: 116.](#)
21. Xia, J. *et al.* (2014) Virus-specific immune memory at peripheral sites of herpes simplex virus type 2 (HSV-2) infection in guinea pigs. [PLoS One. 9 \(12\): e114652.](#)
22. Jeevan A *et al.* (2013) Guinea pig skin, a model for epidermal cellular and molecular changes induced by UVR *in vivo* and *in vitro*: effects on *Mycobacterium bovis* Bacillus Calmette-Guérin vaccination. [Photochem Photobiol. 89 \(1\): 189-98.](#)
23. Gupta A *et al.* (2012) Protective efficacy of *Mycobacterium indicus pranii* against tuberculosis and underlying local lung immune responses in guinea pig model. [Vaccine. 30 \(43\): 6198-209.](#)
24. Gupta, A. *et al.* (2012) Efficacy of *Mycobacterium indicus pranii* immunotherapy as an adjunct to chemotherapy for tuberculosis and underlying immune responses in the lung. [PLoS One. 7 \(7\): e39215.](#)
25. Wu, W.H. *et al.* (2012) Immune status and the development of *Listeria monocytogenes* infection in aged and young guinea pigs. [Clin Invest Med. 35 \(5\): E309.](#)
26. Eckhardt, E. *et al.* (2023) Phosphatidylinositolmannoside vaccination induces lipid-specific Th1-responses and partially protects guinea pigs from *Mycobacterium tuberculosis* challenge. [Sci Rep. 13 \(1\): 18613.](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be

protected from light.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA752F
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA928F\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M445611:250909'

Printed on 09 Sep 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)