

## Datasheet: MCA808F

**BATCH NUMBER 167320**

<b>Description:</b>	MOUSE ANTI RABBIT CD45:FITC
<b>Specificity:</b>	CD45
<b>Other names:</b>	LCA
<b>Format:</b>	FITC
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	L12/201
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rabbit						
<b>Product Form</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid						
<b>Max Ex/Em</b>	<table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation Max (nm)</th> <th>Emission Max (nm)</th> </tr> </thead> <tbody> <tr> <td>FITC</td> <td>490</td> <td>525</td> </tr> </tbody> </table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	FITC	490	525
Fluorophore	Excitation Max (nm)	Emission Max (nm)					
FITC	490	525					
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant						
<b>Buffer Solution</b>	Phosphate buffered saline						
<b>Preservative</b>	0.09% Sodium Azide (NaN <sub>3</sub> )						
<b>Stabilisers</b>	1% Bovine Serum Albumin						
<b>Approx. Protein</b>	IgG concentration 0.1 mg/ml						

## Concentrations

**Immunogen** Glycoproteins isolated from the T cell line, RL-5.

**Fusion Partners** Spleen cells from immunised mice were fused with cells of the P3.X63.Ag8-U1 mouse myeloma cell line.

**Specificity** **Mouse anti Rabbit CD45 antibody, clone L12/201** recognizes the CD45 antigen, also known as leukocyte common antigen (LCA) or T200. Mouse anti Rabbit CD45 antibody, clone L12/201 shows pan leucocyte reactivity by flow cytometry and immunohistochemistry.

Immunoprecipitation was achieved by cross linking antibody to the labelled cell surface yielding a protein migrating by gel electrophoresis at a molecular mass of ~200 kDa.

**Flow Cytometry** Use 10ul of the suggested working dilution to label  $1 \times 10^6$  cells in 100ul

## References

1. Jackson, S. *et al.* (1983) Differentiation antigens identify subpopulations of rabbit T and B lymphocytes. Definition by flow cytometry. [J Exp Med. 157 \(1\): 34-46.](#)
2. Wilkinson, J.M. *et al.* (1984) Cell surface glycoproteins of rabbit lymphocytes: characterization with monoclonal antibodies. [Mol Immunol. 21 \(1\): 95-103.](#)
3. Wilkinson, J.M. *et al.* (1992) A cytotoxic rabbit T-cell line infected with a gamma-herpes virus which expresses CD8 and class II antigens. [Immunology. 77 \(1\): 106-8.](#)
4. Wilkinson, J.M. *et al.* (1993) Immunohistochemical identification of leucocyte populations in normal tissue and inflamed synovium of the rabbit. [J Pathol. 170 \(3\): 315-20.](#)
5. Fenton, M. *et al.* (2001) Cellular senescence after single and repeated balloon catheter denudations of rabbit carotid arteries. [Arterioscler Thromb Vasc Biol. 21: 220-6.](#)
6. Mackenzie, S.M. *et al.* (2006) Immunocontraceptive effects on female rabbits infected with recombinant myxoma virus expressing rabbit ZP2 or ZP3. [Biol Reprod. 74: 511-21.](#)
7. Davis, W.C. & Hamilton, M.J. (2008) Use of flow cytometry to develop and characterize a set of monoclonal antibodies specific for rabbit leukocyte differentiation molecules. [J Vet Sci. 9 \(1\): 51-66.](#)
8. Liang, H. *et al.* (2009) Comparison of the ocular tolerability of a latanoprost cationic emulsion versus conventional formulations of prostaglandins: an *in vivo* toxicity assay. [Mol Vis. 15: 1690-9.](#)
9. Xu, Y. *et al.* (2010) Adenovirus-mediated overexpression of glutathione-s-transferase mitigates transplant arteriosclerosis in rabbit carotid allografts. [Transplantation. 89: 409-16.](#)
10. Vašíček, J *et al.* (2014) Basic Blood Analysis of Rabbits Immunized with Vaccine Against Myxomatosis. [Proc Int Symp Anim Sci 2014: 411-6.](#)
11. Vasicek, J. *et al.* (2015) Determination of Lymphocyte Subset Distribution in the Peripheral Blood of Rabbits Immunized with CFA. [Int Symp Anim Sci UDC:639.112 pp. 226-31.](#)
12. Sijnave, D. *et al.* (2015) Inhibition of Rho-Associated Kinase Prevents Pathological Wound Healing and Neovascularization After Corneal Trauma. [Cornea. 34 \(9\): 1120-9.](#)
13. Ondruska, L. *et al.* (2016) Decrease in C-reactive protein levels in rabbits after vaccination with a live attenuated myxoma virus vaccine [Veterinární Medicína. 61 \(No.](#)

[10\): 571-6.](#)

14. Ondruska, L. *et al.* (2016) Decrease in C-reactive protein levels in rabbits after vaccination with a live attenuated myxoma virus vaccine [Veterinární Medicína. 61 \(No. 10\): 571-6.](#)

15. Kováč, M. *et al.* (2017) Phenotype and ultrastructure of stem cells derived from amniotic fluid of Nitra rabbit [J Cent Eur Agric.18 \(1\): 226-34.](#)

16. Kovac, M. *et al.* (2017) Different RNA and protein expression of surface markers in rabbit amniotic fluid-derived mesenchymal stem cells. [Biotechnol Prog. 33 \(6\): 1601-13.](#)

17. Vašiček, J *et al.* (2018) The Efficiency of Immunomagnetic Sorting of Rabbit Bone Marrow Cell for the Establishment of Mesenchymal Stem Cell Culture. [J Microbiol Biotech Food Sci. 8 \(3\): 890-2.](#)

18. Barth, H. *et al.* (2019) Inflammatory responses after vitrectomy with vitreous substitutes in a rabbit model. [Graefes Arch Clin Exp Ophthalmol. 257 \(4\): 769-83.](#)

19. Kulikova, B. *et al.* (2019) Survivability of rabbit amniotic fluid-derived mesenchymal stem cells post slow-freezing or vitrification. [Acta Histochem. 121 \(4\): 491-9.](#)

20. Desando, G. *et al.* (2021) Prospects on the Potential In Vitro Regenerative Features of Mechanically Treated-Adipose Tissue for Osteoarthritis Care. [Stem Cell Rev Rep. 17 \(4\): 1362-73.](#)

21. Adly, H. *et al.* (2021) Repopulation of multipotent stem cells derived from adult male rabbits on a polycaprolactone scaffold: an *in vitro*. study [Al-Azhar International Medical Journal. 2 \(11\): 37-42.](#)

22. Wee, J.H. *et al.* (2022) Stem cell laden nano and micro collagen/PLGA bimodal fibrous patches for myocardial regeneration. [Biomater Res. 26 \(1\): 79.](#)

23. Mastrolia, I. *et al.* (2022) Autologous Marrow Mesenchymal Stem Cell Driving Bone Regeneration in a Rabbit Model of Femoral Head Osteonecrosis. [Pharmaceutics. 14\(10\):2127.](#)

24. Alpargu, B.T. *et al.* (2025) Effect of TGF- $\beta$ 1 transfected dental pulp derived mesenchymal stem cells on *in vivo* cartilage regeneration. [Iran J Basic Med Sci. 28 \(11\): 1505-1515.](#)

---

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA808F>

---

**Regulatory**

For research purposes only

---

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M430053:240502'

