

Datasheet: MCA911

BATCH NUMBER 158173

Description:	MOUSE ANTI HUMAN CD47
Specificity:	CD47
Other names:	SIRP LIGAND
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	BRIC126
Isotype:	IgG2b
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/10
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting (1)	▪			
Functional Assays (2)	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) **BRIC126 recognizes CD47 under non-reducing conditions, the CD47 epitope recognized is lost on reduction, see [Shahein et al.](#) for details**

(2) **This product contains sodium azide, removal by dialysis is recommended prior to use in functional assays.**

Target Species	Human
Species Cross Reactivity	Reacts with: Dog, Bovine, Sheep, Pig N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for

further information.

Product Form	Purified IgG - liquid
---------------------	-----------------------

Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
--------------------	---

Buffer Solution	TRIS buffered glycine
------------------------	-----------------------

Preservative Stabilisers	0.09% Sodium Azide
---------------------------------	--------------------

Approx. Protein Concentrations	IgG concentration 1 mg/ml
---------------------------------------	---------------------------

Immunogen	Human erythrocytes
------------------	--------------------

External Database Links	UniProt: Q08722 Related reagents Entrez Gene: 961 CD47 Related reagents
--------------------------------	--

Synonyms	MER6
-----------------	------

RRID	AB_321436
-------------	-----------

Specificity	Mouse anti human CD47 antibody, clone BRIC126 recognises the human CD47 cell surface glycoprotein, a heavily N-glycosylated 47-52 kDa molecule. CD47 is expressed on all cells and tissues so far examined, although expression is reduced on erythrocytes of the rare Rh null phenotype.
--------------------	--

Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
-----------------------	---

References	<ol style="list-style-type: none">1. Avent, N. <i>et al.</i> (1988) Monoclonal antibodies that recognize different membrane proteins that are deficient in Rhnull human erythrocytes. One group of antibodies reacts with a variety of cells and tissues whereas the other group is erythroid-specific. Biochem J. 251 (2): 499-505.2. Ide, K. <i>et al.</i> (2007) Role for CD47-SIRPalpha signaling in xenograft rejection by macrophages. Proc Natl Acad Sci U S A. 104 (12): 5062-6.3. Lecchi, C. <i>et al.</i> (2008) Bovine alpha-1 acid glycoprotein can reduce the chemotaxis of bovine monocytes and modulate CD18 expression. Vet Res. 39: 50.4. Chao, M.P. <i>et al.</i> (2011) Therapeutic antibody targeting of CD47 eliminates human acute lymphoblastic leukemia. Cancer Res. 71 (4): 1374-84.5. Shahein, Y.E. <i>et al.</i> (2002) Molecular cloning and functional characterization of the pig homologue of integrin-associated protein (IAP/CD47). Immunology. 106: 564-76.6. Brooke, G. <i>et al.</i> (2004) Human lymphocytes interact directly with CD47 through a novel member of the signal regulatory protein (SIRP) family. J Immunol. 173: 2562-70.
-------------------	--

7. Wewer, C. *et al.* (2011) Transcellular migration of neutrophil granulocytes through the blood-cerebrospinal fluid barrier after infection with *Streptococcus suis*. [J Neuroinflammation. 8: 51.](#)
8. Tichioni, M. *et al.* (2001) Integrin-associated protein (CD47/IAP) contributes to T cell arrest on inflammatory vascular endothelium under flow. [FASEB J. 15: 341-50.](#)
9. Chao, M.P. *et al.* (2010) Anti-CD47 antibody synergizes with rituximab to promote phagocytosis and eradicate non-Hodgkin lymphoma. [Cell. 142: 699-713.](#)
10. Siddhartha, J *et al.* (2016) Methods of manipulating phagocytosis mediated by CD47. [US patent app 15/054930](#)
11. Petrova, P.S. *et al.* (2017) TTI-621 (SIRPαFc): A CD47-Blocking Innate Immune Checkpoint Inhibitor with Broad Antitumor Activity and Minimal Erythrocyte Binding. [Clin Cancer Res. 23 \(4\): 1068-79.](#)

Further Reading	1. Piriou-Guzylack, L. (2008) Membrane markers of the immune cells in swine: an update. Vet Res. 39: 54.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10072 available at: https://www.bio-rad-antibodies.com/SDS/MCA911 10072
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL \(MCA691\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M376223:210123'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)