

## Datasheet: MCA925G

<b>Description:</b>	RAT ANTI MOUSE CD48
<b>Specificity:</b>	CD48
<b>Other names:</b>	BLAST-1
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	OX-78
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.25 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide

<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1 mg/ml
<b>Immunogen</b>	Recombinant soluble mouse CD48.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P18181</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">12506</a> Cd48    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	Bcm-1
<b>RRID</b>	AB_321439
<b>Fusion Partners</b>	Spleen cells from immunized AO Rats were fused with cells from the Rat Y3 myeloma cell line.
<b>Specificity</b>	<p><b>Rat anti Mouse CD48 antibody, clone OX-78</b> recognizes murine CD48 antigen, also known as BCM1 surface antigen, BLAST-1, HM48-1, MRC OX-45 surface antigen, SLAM family member 2. CD48 is a 195 amino acid GPI anchored membrane glycoprotein with an additional 22 amino acid signal peptide and a 23 amino acid propeptide region cleaved in the mature form.</p> <p>Rat anti Mouse CD48 antibody, clone OX-78 labels mouse thymocytes, B and T cells and macrophages. CD48 is a ligand for the mouse CD2 cell surface antigen.</p>
<b>Flow Cytometry</b>	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells in 100ul.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Kato, K. <i>et al.</i> (1992) CD48 is a counter-receptor for mouse CD2 and is involved in T cell activation. <a href="#">J Exp Med. 176 (5): 1241-9.</a></li> <li>2. Lapaque, N. <i>et al.</i> (2006) Characterization of <i>Brucella abortus</i> lipopolysaccharide macrodomains as mega rafts. <a href="#">Cell Microbiol. 8 (2): 197-206.</a></li> <li>3. Kashiwakura, Y. <i>et al.</i> (2013) CD2-mediated regulation of peripheral CD4(+) CD25(+) regulatory T-cell apoptosis accompanied by down-regulation of Bim. <a href="#">Immunology. 139 (1): 48-60.</a></li> <li>4. Wild, M.K. <i>et al.</i> (1999) Dependence of T cell antigen recognition on the dimensions of an accessory receptor-ligand complex. <a href="#">J Exp Med. 190: 31-41.</a></li> <li>5. Grassinger, J. <i>et al.</i> (2010) Phenotypically identical hemopoietic stem cells isolated from different regions of bone marrow have different biologic potential. <a href="#">Blood. 116 (17): 3185-96.</a></li> <li>6. Oh, J. <i>et al.</i> (2018) MARCH1 protects the lipid raft and tetraspanin web from MHCII proteotoxicity in dendritic cells. <a href="#">J Cell Biol. 217 (4): 1395-410.</a></li> </ol>
<b>Storage</b>	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for

short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: 10040: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a>
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rat IgG (STAR69...)	<a href="#">FITC</a>
Goat Anti Rat IgG (STAR73...)	<a href="#">RPE</a>
Rabbit Anti Rat IgG (STAR17...)	<a href="#">FITC</a>
Goat Anti Rat IgG (STAR131...)	<a href="#">Alk. Phos.</a> , <a href="#">Biotin</a>
Goat Anti Rat IgG (STAR72...)	<a href="#">HRP</a>
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...)	<a href="#">DyLight®550</a> , <a href="#">DyLight®650</a> , <a href="#">DyLight®800</a>
Rabbit Anti Rat IgG (STAR21...)	<a href="#">HRP</a>
Rabbit Anti Rat IgG (STAR16...)	<a href="#">DyLight®800</a>

### Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL \(MCA1212\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M381905:210512'

Printed on 10 Sep 2022