

Datasheet: MPP017

BATCH NUMBER 156268

Description:	NEISSERIA GONORRHOEAE
Name:	NEISSERIA GONORRHOEAE
Format:	Inactivated Pathogen
Product Type:	Antigen
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Inactivated, homogenised <i>Neisseria gonorrhoeae</i> - frozen
Preparation	A homogenised pure bacteria. The preparation is inactivated by heat treatment at 56°C for 3 hours. Inactivation is confirmed by no visible colonies on permissive plates after 5 days.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Total protein concentration 0.37 mg/ml

Product Information *Neisseria gonorrhoeae* is a non-motile, Gram-negative aerobic coccus with a typical outer membrane composed of proteins, phospholipids and lipopolysaccharide (LPS). Neisserial lipopolysaccharide has a highly branched basal oligosaccharide structure and no

repeating O-antigen subunits, also referred to as lipooligosaccharide (LOS).

Neisseria gonorrhoeae is responsible for the sexually transmitted infection gonorrhoea, a type of urethritis. The infection is characterised by a purulent discharge, but asymptomatic infection is a major problem. *Neisseria gonorrhoeae* could also cause conjunctivitis, pharyngitis or proctitis. The infection can be treated with antibiotic effective against resistant strains.

Activity	This product has been rendered inactive by standard procedures. However this material should still be handled as infectious and should be disposed of appropriately.
-----------------	--

Storage	Store at -70°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10258 available at: https://www.bio-rad-antibodies.com/SDS/MPP017 10258
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M369242:200529'

Printed on 27 Apr 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)