

Datasheet: MPP024

BATCH NUMBER 167509

Description:	CLOSTRIDIUM DIFFICILE TOXOID B
Name:	CLOSTRIDIUM DIFFICILE TOXOID B
Format:	Purified
Product Type:	Antigen
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Inactivated toxoid derived from purified <i>Clostridium difficile</i> Toxin B - lyophilized
Reconstitution	Reconstitute with 250 µl distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Inactivated using formaldehyde. No cell rounding in the sensitive verocytotoxicity assay at 2 µg/ml compared to active toxin lethal dose of 4 pg/ml.
Buffer Solution	0.05M HEPES, 0.15M NaCl, 5% Sucrose
Preservative Stabilisers	None present
Approx. Protein Concentrations	0.2 mg/ml after reconstitution
Product Information	<i>Clostridium difficile</i> is a gram-positive bacteria naturally present in the gut of some

children and very few adults. Occasionally, the use of antibiotics to treat other infections may lead to *C. difficile* infection and severe diarrhoea and fever. Symptoms are caused by toxin A (an enterotoxin) and toxin B (a cytotoxin).

MPP024 is a highly purified, non toxic preparation derived from *Clostridium difficile* Toxin B. Antigenicity of this preparation is maintained to be no less than 90% of original toxin strength as measured by direct ELISA.

For a preparation of Toxoid A, see [MPP022](#).

Purity	>95% by SDS-PAGE
Storage	<p>This product is shipped at ambient temperature. Prior to reconstitution store at +4°C. After reconstitution it is recommended to aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the product. Storage in frost-free freezers is not recommended.</p>
Guarantee	1 month from date of reconstitution if stored at +4°C
Health And Safety Information	Material Safety Datasheet documentation #10505 available at: https://www.bio-rad-antibodies.com/SDS/MPP024 10505
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M418660:230427'

Printed on 18 Jan 2024