

Datasheet: OBT0943P

BATCH NUMBER 149794

Description:	RABBIT ANTI LEGIONELLA SPP:HRP
Specificity:	LEGIONELLA SPP
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			1/10 - 1/100
Immunohistology - Paraffin			▪	
ELISA	▪			1/100 - 1/500
Immunoprecipitation			▪	
Immunoblotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Ig Fraction conjugated to Horseradish Peroxidase (HRP) - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.002% Thiomersal 1.0% w/v Bovine Serum Albumin
Approx. Protein Concentrations	1 mg/ml
Immunogen	Whole cell preparation of <i>Legionella pneumophila</i> ATCC#33152

RRID AB_609812

Specificity Rabbit anti **Legionella SPP antibody** recognizes *Legionella* spp, recognizing multiple antigens, and has been shown to be reactive with serogroups 1-12.

Rabbit anti Legionella SPP antibody is not absorbed, so may cross-react with related microorganisms.

References

1. Yáñez, M.A. *et al.* (2005) Quantitative detection of Legionella pneumophila in water samples by immunomagnetic purification and real-time PCR amplification of the dotA gene. [Appl Environ Microbiol. 71: 3433-41.](#)
2. Allegra, S. *et al.* (2011) Evaluation of an immunomagnetic separation assay in combination with cultivation to improve Legionella pneumophila serogroup 1 recovery from environmental samples. [J Appl Microbiol. 110: 952-961.](#)

Storage Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10095 available at: <https://www.bio-rad-antibodies.com/SDS/OBT0943P10095>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M369283:200529'

Printed on 18 Jan 2024

