

Datasheet: OBT0947

BATCH NUMBER 165071

Description:	RABBIT ANTI MYCOBACTERIUM TUBERCULOSIS
Specificity:	MYCOBACTERIUM TUBERCULOSIS
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Purified IgG - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 5.0 mg/ml
Immunogen	Purified protein derivative (PPD) from <i>M. tuberculosis</i>
RRID	AB_620694

Specificity Rabbit anti *Mycobacterium tuberculosis* polyclonal antibody recognizes PPD from *Mycobacterium tuberculosis*. Rabbit anti *M. tuberculosis* has not been cross absorbed and may react with related micro-organisms, however the antibody is non-reactive with *E.coli* K12, *Salmonella typhimurium*, *Pseudomonas aeruginosa*, *Streptococcus* (group B), *Candida albicans* and *Neisseria meningitidis*.

- References**
1. Ferrer, N.L. *et al.* (2009) Intracellular replication of attenuated *Mycobacterium tuberculosis* phoP mutant in the absence of host cell cytotoxicity. [Microbes Infect. 11: 115-22.](#)
 2. Sanchez, J. *et al.* (2011) Microscopical and Immunological Features of Tuberculoid Granulomata and Cavitory Pulmonary Tuberculosis in Naturally Infected Goats. [J Comp Pathol. 145: 107-17.](#)
 3. Buendía, A.J. *et al.* (2013) Ante-mortem diagnosis of caprine tuberculosis in persistently infected herds: influence of lesion type on the sensitivity of diagnostic tests. [Res Vet Sci. 95 \(3\): 1107-13.](#)
 4. Herrtwich, L. *et al.* (2016) DNA Damage Signaling Instructs Polyploid Macrophage Fate in Granulomas. [Cell. 167 \(5\): 1264-1280.e18.](#)
 5. BernardinSouibgui, C. *et al.* (2017) Virulence test using nematodes to prescreen *Nocardia* species capable of inducing neurodegeneration and behavioral disorders. [PeerJ. 5: e3823.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/OBT0947>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

'M414080:221130'

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)