

Datasheet: OBT1424

Description:	SHEEP ANTI HUMAN 17-ALPHA-HYDROXYPROGESTERONE-3-CMO
Specificity:	17-ALPHA-HYDROXYPROGESTERONE-3-CMO
Other names:	17OHP
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Serum - liquid
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
Immunogen	17-Alpha-Hydroxyprogesterone-3-CMO conjugated to KLH.
RRID	AB_618537

Specificity

Sheep anti Human 17-alpha-Hydroxyprogesterone-3-Cmo antibody recognizes 17-Alpha-Hydroxyprogesterone-3-CMO, a progestogen which may reduce the risk of pre-term delivery in pregnant women.

This product has the following reactivities:

17-Alpha-Hydroxyprogesterone	100%
11-Deoxycortisol	<2%
Progesterone	<1%

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10081 available at:
10081: <https://www.bio-rad-antibodies.com/uploads/MSDS/10081.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [DyLight®488](#), [HRP](#)

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M389755:210806'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)