

## Datasheet: OBT1424

<b>Description:</b>	SHEEP ANTI HUMAN 17-ALPHA-HYDROXYPROGESTERONE-3-CMO
<b>Specificity:</b>	17-ALPHA-HYDROXYPROGESTERONE-3-CMO
<b>Other names:</b>	17OHP
<b>Format:</b>	Serum
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Serum - liquid
<b>Preservative Stabilisers</b>	<0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Immunogen</b>	17-Alpha-Hydroxyprogesterone-3-CMO conjugated to KLH.
<b>RRID</b>	AB_618537

**Specificity** **Sheep anti Human 17-alpha-Hydroxyprogesterone-3-Cmo antibody** recognizes 17-Alpha-Hydroxyprogesterone-3-CMO, a progestogen which may reduce the risk of pre-term delivery in pregnant women.

This product has the following reactivities:

17-Alpha-Hydroxyprogesterone	100%
11-Deoxycortisol	<2%
Progesterone	<1%

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

**Guarantee** 12 months from date of despatch

**Health And Safety Information** Material Safety Datasheet documentation #10081 available at: <https://www.bio-rad-antibodies.com/SDS/OBT1424>  
10081

**Regulatory** For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376	<b>Worldwide</b>	Tel: +44 (0)1865 852 700	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21
	Fax: +1 919 878 3751		Fax: +44 (0)1865 852 739		Fax: +49 (0) 89 8090 95 50
	Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>		Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>		Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M429080:240315'

Printed on 15 Mar 2024