

Datasheet: OBT1521F

BATCH NUMBER 170515

Description:	RABBIT ANTI MOUSE IgG2a:FITC
Specificity:	IgG2a
Format:	FITC
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse						
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - lyophilised						
Reconstitution	Reconstitute with 1.0ml distilled water						
Max Ex/Em	<table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation Max (nm)</th> <th>Emission Max (nm)</th> </tr> </thead> <tbody> <tr> <td>FITC</td> <td>490</td> <td>525</td> </tr> </tbody> </table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	FITC	490	525
Fluorophore	Excitation Max (nm)	Emission Max (nm)					
FITC	490	525					
Preparation	Purified IgG prepared by affinity chromatography.						
Buffer Solution	Phosphate buffered saline						
Preservative	0.01% Sodium azide						
Stabilisers	1% w/v Bovine serum albumin						
Approx. Protein Concentrations	1.0 mg/ml						

Immunogen Mouse IgG2a heavy chain

External Database

Links

UniProt:

[P01865](#) [Related reagents](#)

[P01863](#) [Related reagents](#)

[P01864](#) [Related reagents](#)

Entrez Gene:

[380793](#) Igh-1a [Related reagents](#)

[380793](#) Igh-1a [Related reagents](#)

[380793](#) Igh-1a [Related reagents](#)

RRID AB_616783

Specificity **Rabbit anti Mouse IgG2a antibody** recognises murine IgG2a and shows <1% cross-reactivity with other murine immunoglobulin classes or subclasses. It also shows <1% cross-reactivity with human heavy and light chains.

References 1. Paillot, R. *et al.* (2010) Contribution of each of four Superantigens to *Streptococcus equi*-induced mitogenicity, gamma interferon synthesis, and immunity. [Infect Immun. 78: 1728-39.](#)

Storage Prior to reconstitution store at +4°C.
After reconstitution store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/OBT1521F>
10041

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M369320:200529'

Printed on 18 Jan 2024
