

Datasheet: OBT1646

BATCH NUMBER 163439

| | |
|----------------------|--------------------------------------|
| Description: | CHICKEN ANTI INFLUENZA HEMAGGLUTININ |
| Specificity: | INFLUENZA HEMAGGLUTININ |
| Other names: | HA |
| Format: | Purified |
| Product Type: | Polyclonal Antibody |
| Isotype: | IgY |
| Quantity: | 0.1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | ▪ | | | 1/250 |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | | | ▪ | |
| Immunoprecipitation | ▪ | | | 1/250 |
| Western Blotting | ▪ | | | 1/500 |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|--|
| Target Species | Viral |
| Product Form | Purified IgG - liquid |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 0.02% Sodium Azide (NaN ₃) |
| Approx. Protein Concentrations | IgG concentration 0.1 mg/ml |
| Immunogen | Synthetic peptide YPYDVDPYA from influenza hemagglutinin (HA). |

| | |
|--------------------------------------|---|
| RRID | AB_620520 |
| Specificity | Chicken anti influenza heamagglutinin antibody recognizes influenza hemagglutinin, a glycosylated antireceptor which binds to cell surface glycoproteins prior to penetration of the cell membrane and uncoating of the virus. |
| Storage | Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use. |
| Guarantee | 6 months from date of despatch |
| Health And Safety Information | Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/OBT1646 10040 |
| Regulatory | For research purposes only |

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M374632:201120'

Printed on 23 May 2025