

Datasheet: OBT1680

Description:	MOUSE ANTI HUMAN CYSTIC FIBROSIS TRANSMEMBRANE CONDUCTANCE REGULATOR
Specificity:	CFTR
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	M3A7
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen (1)	▪			
Immunohistology - Paraffin		▪		
Immunoprecipitation	▪			
Western Blotting (2)	▪			10ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)The epitope recognised by this antibody is reported to be sensitive to formaldehyde fixation and tissue processing. Bio-Rad recommends the use of acetone fixation for frozen sections.

(2)This antibody recognises a 170kDa band under native conditions.

Target Species	Human
Species Cross Reactivity	Does not react with:Mouse
Product Form	Purified IgG - liquid
Preparation	Antibody purified from ascites
Buffer Solution	Phosphate Buffered Saline

Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Fusion protein containing residues 1197-1480 of human CFTR (Chang et al. 1993).
External Database Links	<p>UniProt: P13569 Related reagents</p> <p>Entrez Gene: 1080 CFTR Related reagents</p>
Synonyms	ABCC7
RRID	AB_617226
Specificity	Mouse anti Human Cystic Fibrosis Transmembrane Conductance Regulator antibody, clone M3A7 recognizes the human cystic fibrosis transmembrane conductance regulator (CFTR). It recognises an epitope localizing to the ABC transporter 2 domain, specifically to a region between amino acids 1370 to 1380.
Protein Molecular Weight	170kDA
References	<ol style="list-style-type: none"> 1. Kartner, N. <i>et al.</i> (1992) Mislocalization of delta F508 CFTR in cystic fibrosis sweat gland. Nat Genet. 1 (5): 321-7. 2. Chang, X.B. <i>et al.</i> (1993) Protein kinase A (PKA) still activates CFTR chloride channel after mutagenesis of all 10 PKA consensus phosphorylation sites. J Biol Chem. 268 (15): 11304-11. 3. Kartner, N. & Riordan, J.R. (1998) Characterization of polyclonal and monoclonal antibodies to cystic fibrosis transmembrane conductance regulator. Methods Enzymol. 292: 629-52. 4. Zhang, F. <i>et al.</i> (1998) Limited proteolysis as a probe for arrested conformational maturation of delta F508 CFTR. Nat Struct Biol. 5 (3): 180-3. 5. Benharouga, M. <i>et al.</i> (2001) COOH-terminal truncations promote proteasome-dependent degradation of mature cystic fibrosis transmembrane conductance regulator from post-Golgi compartments. J Cell Biol. 153: 957-70.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/OBT1680>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...) [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M406014:220916'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)