

Datasheet: OBT1753

Description:	RABBIT ANTI HUMAN DRP-3
Specificity:	DRP-3
Other names:	DPYSL3
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunoprecipitation	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Reacts with: Rat, Cat, Mouse, Monkey N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Serum - liquid
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Immunogen	Peptide with the sequence YDGPVFDLTTTPK.
External Database Links	<p>UniProt: Q14195 Related reagents</p> <p>Entrez Gene: 1809 DPYSL3 Related reagents</p>
Synonyms	CRMP4, DRP3, ULIP, ULIP1
RRID	AB_620371

Specificity **Rabbit anti Human DRP-3 antibody** recognizes human Dihydropyrimidinase-related protein 3, also known as DRP-3, Collapsin response mediator protein 4, TUC-4, CRMP-4, Unc-33-like phosphoprotein 1 or Ulip-1. DRP-3 is a 570 amino acid ~60 kDa cytoplasmic protein expressed in neuronal cells. Rabbit anti Human DRP-3 antibody does not recognise TUC-1, TUC-2, or TUC-3 proteins.

Further Reading 1. Moreira, J.M. *et al.* (2010) Tissue proteomics of the human mammary gland: towards an abridged definition of the molecular phenotypes underlying epithelial normalcy. [Mol Oncol. 4 \(6\): 539-61.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10342 available at: 10342: <https://www.bio-rad-antibodies.com/uploads/MSDS/10342.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@680](#), [DyLight@800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M369349:200529'

Printed on 11 Aug 2020

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)