

Datasheet: OBT2003

BATCH NUMBER 156151

Description:	MOUSE ANTI KINETOPLASTID MEMBRANE PROTEIN 11
Specificity:	KINETOPLASTID MEMBRANE PROTEIN 11
Other names:	KMP-11
Format:	Ascites
Product Type:	Monoclonal Antibody
Clone:	L-157
Isotype:	IgG1
Quantity:	0.5 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/2000
Immunofluorescence	▪			1/500 - 1/1000
Immunoblotting	▪			1/1000 - 1/3000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Protozoan
Product Form	Ascites - lyophilised
Reconstitution	Reconstitute with 0.5 ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preservative Stabilisers	None present
Immunogen	<i>L. donovani</i> lipophosphoglycan semi purified from epimastigotes.

**External Database
Links**

UniProt:

[Q36736](#)

[Related reagents](#)

RRID

AB_619092

Fusion Partners

Spleen cells from immunised BALB/c mice were fused with cells of the SP2/0 myeloma cell line.

Specificity

Mouse anti Kinetoplastid Membrane Protein 11 antibody, clone L-157 recognizes Kinetoplastid Membrane Protein-11 (KMP-11), previously known as the lipophosphoglycan-associated protein (LPGAP), a 92 amino acid ~11 kDa membrane protein. Mouse anti Kinetoplastid Membrane Protein 11 antibody, clone L-157 was originally derived against *Leishmania donovani* lipophosphoglycan, but found to recognize the tightly associated LPGAP. KMP-11 is found in a wide variety of Kinetoplastid parasites including African trypanosomes, *Leishmania* and *Crithidia* ([Tolson et al. 1994](#)).

References

1. Tolson, D.L. *et al.* (1994) The kinetoplastid membrane protein 11 of *Leishmania donovani* and African trypanosomes is a potent stimulator of T-lymphocyte proliferation. [Infect Immun. 62 \(11\): 4893-9.](#)
2. Pearson, T.W. *et al.* (2000) The major cell surface glycoprotein procyclin is a receptor for induction of a novel form of cell death in African trypanosomes *in vitro*. [Mol Biochem Parasitol. 111 \(2\): 333-49.](#)
3. Jardim, A. *et al.* (2002) Peroxisomal targeting protein 14 (PEX14) from *Leishmania donovani*. Molecular, biochemical, and immunocytochemical characterization. [Mol Biochem Parasitol. 124 \(1-2\): 51-62.](#)
4. McNeely, T.B. *et al.* (1990) Characterization of *Leishmania donovani* variant clones using anti-lipophosphoglycan monoclonal antibodies. [Glycobiology. 1 \(1\): 63-9.](#)
5. Jardim, A. *et al.* (1991) The *Leishmania donovani* lipophosphoglycan T lymphocyte-reactive component is a tightly associated protein complex. [J Immunol. 147 \(10\): 3538-44.](#)

Storage

Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10484 available at:
10484: <https://www.bio-rad-antibodies.com/uploads/MSDS/10484.pdf>

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

Printed on 08 Feb 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)