

Datasheet: PAP004

BATCH NUMBER 159058

Description:	RECOMBINANT CHICKEN INTERFERON ALPHA
Name:	IFN ALPHA
Other names:	INTERFERON ALPHA
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA			▪	
Functional Assays	▪			

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Chicken
Product Form	Purified recombinant protein expressed in <i>E.coli</i> - liquid.
Preparation	Recombinant protein prepared from <i>E.coli</i>
Buffer Solution	0.05M Sodium Acetate 1M Sodium Chloride 0.01M Tris-HCl 2M Urea 10mM Beta-Mercaptoethanol
Preservative Stabilisers	None present
Approx. Protein Concentrations	Total protein concentration 1.0 mg/ml

**External Database
Links**

UniProt:

[P42165](#) [Related reagents](#)

Entrez Gene:

[396398](#) IFNA3 [Related reagents](#)

Product Information **Recombinant Chicken interferon alpha** (IFN alpha) shares structural and biological properties with mammalian IFN1. Recombinant chicken IFN alpha (rChIFNa) has been shown to act as a highly potent anti-viral agent *in vitro* ([Schultz, U. et al. 1995](#)). rChIFNa also resulted in a noticeable delay in the progression of tumors in CC progressor chickens infected with Rous Sarcoma Virus, when administered intravenously ([Plachy, J. et al. 1999](#)).

Recombinant chicken interferon alpha is the immunogen used to produce the polyclonal Rabbit anti Chicken [interferon alpha antibody](#).

**Protein Molecular
Weight**

19 kD

Activity

1x10⁷units/mg

Purity

>85% by Ni chelate chromatography

References

1. Schultz, U. *et al.* (1995) Recombinant chicken interferon: a potent antiviral agent that lacks intrinsic macrophage activating factor activity. [Eur J Immunol. 25 \(3\): 847-51.](#)
2. Plachy, J. *et al.* (1999) Protective effects of Type I and Type II interferons toward Rous Sarcoma Virus-induced tumors in chickens. [Virology. 256: 85-91.](#)
3. Schwarz, H. *et al.* (2004) Synthesis of IFN-beta by virus-infected chicken embryo cells demonstrated with specific antisera and a new bioassay. [J Interferon Cytokine Res. 24 \(3\): 179-84.](#)
4. Schultz, U. *et al.* (1995) Recombinant chicken interferon from *Escherichia coli* and transfected COS cells is biologically active. [Eur J Biochem. 229 \(1\): 73-6.](#)
5. Jiang, H. *et al.* (2011) Chicken interferon alpha pretreatment reduces virus replication of pandemic H1N1 and H5N9 avian influenza viruses in lung cell cultures from different avian species. [Virology. 8: 447.](#)
6. Smith, S.E. *et al.* (2013) Chicken interferon-inducible transmembrane protein 3 restricts influenza viruses and lyssa viruses *in vitro*. [J Virol. 87 \(23\): 12957-66.](#)
7. Alamares, J.G. *et al.* (2010) The interferon antagonistic activities of the V proteins from two strains of Newcastle disease virus correlate with their known virulence properties. [Virus Res. 147 \(1\): 153-7.](#)
8. Lostalé-seijo, I. *et al.* (2016) Interferon induction by avian reovirus. [Virology. 487: 104-11.](#)
9. Lostalé-seijo, I. *et al.* (2016) Response of Three Different Viruses to Interferon Priming and Dithiothreitol Treatment of Avian Cells. [J Virol. 90 \(18\): 8328-40.](#)
10. Deng, L. *et al.* (2016) Characterization and functional studies of fowl adenovirus 9 dUTPase. [Virology. 497: 251-61.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein.
This product should be stored undiluted.
Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 3 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #20481 available at:
<https://www.bio-rad-antibodies.com/SDS/PAP004>
20481

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373928:201014'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)