

Datasheet: PBP005

BATCH NUMBER 166982

Description:	RECOMBINANT BOVINE TNF ALPHA
Name:	TNF ALPHA
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	25 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Functional Assays	▪			

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only.

Target Species	Bovine
Product Form	Purified Recombinant Protein- lyophilized
Reconstitution	Reconstitute with 0.5 ml distilled water. Further dilutions should be made in a buffer containing carrier protein. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i>
Source	<i>Pichia pastoris</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present

**External Database
Links**

UniProt:

[Q06599](#) [Related reagents](#)

Entrez Gene:

[280943](#) TNF [Related reagents](#)

Synonyms TNFA, TNFSF2

Product Information **Bovine Tumour Necrosis Factor alpha (TNF alpha)** is produced by a range of cells, including macrophages, monocytes, T and B lymphocytes. TNF alpha is an important mediator of the bovine inflammatory response.

Bovine Tumour Necrosis Factor alpha (**PBP005**) has an ED₅₀ in the range of 10-25 ng/ml, determined by measuring the cytolytic effect of recombinant bovine TNF alpha on murine WEHI-13VAR cells in the presence of Actinomycin D.

Protein Molecular Weight 17.4kD

Purity >95% by SDS PAGE analysis

ELISA PBP005 may be used as a standard with [MCA2334](#) or [MCA2335B](#).

References 1. Simojoki, H. *et al.* (2011) Innate immune response in experimentally induced bovine intramammary infection with *Staphylococcus simulans* and *S. epidermidis*. [Vet Res. 42: 49.](#)

Storage Prior to reconstitution store at +4°C.
Following reconstitution store at -70°C.
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10302 available at:
<https://www.bio-rad-antibodies.com/SDS/PBP005>
10302

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M426021:231101'

Printed on 19 Jan 2024

