

Datasheet: PBP017

BATCH NUMBER 154101

Description:	PURIFIED BOVINE FIBRONECTIN
Name:	FIBRONECTIN
Format:	Purified
Product Type:	Purified Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Tissue Culture	▪			1 - 10 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified bovine fibronectin - sterile liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Endotoxin Level	<10 EU/mg
Approx. Protein Concentrations	Protein concentration 0.5 mg/ml

External Database Links

UniProt:

[P07589](https://www.uniprot.org/entry/P07589) [Related reagents](#)

Entrez Gene:

[280794](https://www.ncbi.nlm.nih.gov/gene/280794) FN1 [Related reagents](#)

Product Information	Fibronectin is a 450 kDa glycoprotein involved in cell adhesion and migration processes including embryogenesis, metastasis, wound healing, blood coagulation and host immunity. Fibronectin may be used as a coating on tissue culture surfaces to promote attachment, spreading and proliferation of various cell types.
Purity	>=95% by SDS PAGE
Instructions For Use	<p>Preparation of Fibronectin:</p> <p>Generally fibronectin should be stored frozen or refrigerated at 1-5 mg/ml. If frozen, thaw at 37°C. Fibronectin should never be vortexed or treated roughly as it precipitates out of solution easily (visible as a granulated precipitate or as a slimy strand or clumps). Any observed insoluble aggregates are very difficult to resuspend. If precipitation is observed, the solution should initially be warmed to 37°C with gentle shaking, followed by aseptic filtration via a 5 µm (or smaller) cut off syringe filter.</p> <p>It is also recommended that fibronectin should not be stored long term in buffer solutions containing Mg²⁺ or Ca²⁺ as these can contribute to fibronectin precipitation over time.</p> <p>General Coating Procedure (for 96-well plates):</p> <ol style="list-style-type: none"> 1. Plates: Choice of plates can affect the amount of protein successfully coated. Bio-Rad recommends the use of commercially available high protein binding plates. 2. Add between 1-10 µg/ml of fibronectin in PBS buffer for an overnight (or longer) incubation at 2-8°C (100 µl/well). The optimum amount of protein to saturate the plate surface should be determined by the user. For ELISA applications some users have coated with carbonate buffer pH 9.0 but this is generally not necessary and can damage the cell binding activity of the fibronectin. 3. Regardless of assay type (cell binding, ELISA etc.), it is necessary to block the remaining protein adherence sites on the plate. Therefore, in a separate step add a blocking protein such as 2-5% BSA in PBS for at least 1 hr at room temperature or overnight at 2-8°C (200 µl/well). Caution: The BSA solution should be 0.2 µm filtered prior to use to remove excess non-specific binding effects in the assay caused by insoluble BSA clumps. The use of RIA-grade (radio-immunoassay grade, suitable for immunoassays) BSA is recommended. 4. Blocked plates can stored in the refrigerator for several weeks or can be decanted and dried and stored for months in a dessicator. Dessicated material should be rehydrated for 15 min with PBS before use.
Storage	<p>Store at +4°C or at -20°C for long term storage.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.</p>
Guarantee	12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10209 available at:
<https://www.bio-rad-antibodies.com/SDS/PBP017>
10209

Regulatory

For research purposes only

**North & South
America** Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M369473:200529'

Printed on 01 May 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)