

Datasheet: PBP025

**BATCH NUMBER 151725**

<b>Description:</b>	RECOMBINANT BOVINE INTERLEUKIN-17A
<b>Name:</b>	IL-17A
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	5 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Bovine
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	10% Trehalose
<b>External Database Links</b>	<b>UniProt:</b> <a href="https://www.uniprot.org/entry/Q687Y7">Q687Y7</a> <a href="#">Related reagents</a>

**Entrez Gene:**[282863](#) IL17A [Related reagents](#)

---

**Product Information** **Recombinant Bovine interleukin-17A** is a recombinant protein corresponding to aa 22 - 153 of bovine interleukin-17A (IL-17A). IL-17A is a disulfide-linked homodimer and member of the IL-17 family, secreted primarily by activated CD4<sup>+</sup> T cells, which signal through the IL-17 receptor (IL-17R).

IL-17 induces a variety of target cells to secrete inflammatory cytokines, including IL-1 $\beta$ , IL-6, TNF $\alpha$ , IFN $\gamma$  and granulocyte colony-stimulating factor, as well as synergizing with other inflammatory cytokines and agonists.

---

**Protein Molecular Weight** 15.0 kDa

---

**Purity** >95% by SDS PAGE analysis

---

**ELISA** This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2376](#) as a coating antibody and [AHP2376B](#) as a detection antibody.

---

**Storage** Prior to reconstitution store at -20°C.  
After reconstitution store at -20°C.  
Storage in frost-free freezers is not recommended. This product should be stored undiluted.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10591 available at: <https://www.bio-rad-antibodies.com/SDS/PBP025>  
10591

---

**Regulatory** For research purposes only

---

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M334890:181203'

**Printed on 19 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)