

Datasheet: PBP032

BATCH NUMBER 158030

| | |
|----------------------|------------------------------------|
| Description: | RECOMBINANT BOVINE INTERFERON BETA |
| Name: | IFN BETA |
| Format: | Rec. Protein |
| Product Type: | Recombinant Protein |
| Quantity: | 5 µg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |
| Western Blotting | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|--|
| Target Species | Bovine |
| Product Form | Purified recombinant protein - lyophilized |
| Reconstitution | Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. |
| Preparation | Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 10% Trehalose |
| Product Information | Recombinant Bovine interferon beta (IFN β) is a recombinant protein of approximately 19.4 kDa consisting of a sequence of 164 aa. |

Bovine IFN β is a type I interferon belonging to alpha/beta interferon family and involved in defense against viral and bacterial infections.

| | |
|---------------------------------|----------|
| Protein Molecular Weight | 19.4 kDa |
|---------------------------------|----------|

| | |
|---------------|---------------------------|
| Purity | >95% by SDS PAGE analysis |
|---------------|---------------------------|

| | |
|--------------|--|
| ELISA | This reagent may be used as a standard in a sandwich ELISA assay in conjunction with AHP2371 as a coating antibody and AHP2371B as a detection antibody. |
|--------------|--|

| | |
|----------------|--|
| Storage | Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. |
|----------------|--|

| | |
|------------------|---------------------------------|
| Guarantee | 12 months from date of despatch |
|------------------|---------------------------------|

| | |
|--------------------------------------|---|
| Health And Safety Information | Material Safety Datasheet documentation #10591 available at: https://www.bio-rad-antibodies.com/SDS/PBP032 10591 |
|--------------------------------------|---|

| | |
|-------------------|----------------------------|
| Regulatory | For research purposes only |
|-------------------|----------------------------|

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M334895:181203'

Printed on 15 Jul 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)