

Datasheet: PBP039

**BATCH NUMBER 161661**

<b>Description:</b>	RECOMBINANT BOVINE INTERLEUKIN-8
<b>Name:</b>	IL-8
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	5 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Bovine
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	10% Trehalose
<b>External Database Links</b>	<b>UniProt:</b> <a href="https://www.uniprot.org/entry/P79255">P79255</a> <a href="#">Related reagents</a>

**Entrez Gene:**[280828](#) IL8 [Related reagents](#)

---

<b>Synonyms</b>	CXCL8
-----------------	-------

---

<b>Product Information</b>	<b>Recombinant Bovine Interleukin-8 (IL-8)</b> is a recombinant protein corresponding to aa 23-101 of native bovine IL-8.  Bovine IL-8 is produced by macrophages and released from several cell types in response to an inflammatory stimulus. IL-8 is a member of the chemokine family and as chemotactic factor attracts neutrophils, basophils, and T cells but not monocytes.
----------------------------	--

---

<b>Protein Molecular Weight</b>	9.1 kDa
---------------------------------	---------

---

<b>Purity</b>	>95% by SDS PAGE analysis
---------------	---------------------------

---

<b>ELISA</b>	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with <a href="#">AHP2817</a> as a coating antibody and <a href="#">AHP2817B</a> as a detection antibody.
--------------	--

---

<b>Storage</b>	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
----------------	--

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10591 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PBP039">https://www.bio-rad-antibodies.com/SDS/PBP039</a> 10591
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373617:200923'

Printed on 19 Jan 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)