

Datasheet: PBP039

Description:	RECOMBINANT BOVINE INTERLEUKIN-8
Name:	IL-8
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose

External Database Links

UniProt:

[P79255](#)

[Related reagents](#)

Entrez Gene:

Synonyms	CXCL8
-----------------	-------

Product Information	Recombinant Bovine Interleukin-8 (IL-8) is a recombinant protein corresponding to aa 23-101 of native bovine IL-8. Bovine IL-8 is produced by macrophages and released from several cell types in response to an inflammatory stimulus. IL-8 is a member of the chemokine family and as chemotactic factor attracts neutrophils, basophils, and T cells but not monocytes.
----------------------------	--

Protein Molecular Weight	9.1 kDa
---------------------------------	---------

Purity	>95% by SDS PAGE analysis
---------------	---------------------------

ELISA	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with AHP2817 as a coating antibody and AHP2817B as a detection antibody.
--------------	--

Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10591 available at: 10591: https://www.bio-rad-antibodies.com/uploads/MSDS/10591.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M373617:200923'

Printed on 21 Mar 2022