

## Datasheet: PDP001

<b>Description:</b>	RECOMBINANT DOG TNF ALPHA
<b>Name:</b>	TNF ALPHA
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	20 µg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA			▪	
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Dog
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 0.5 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.  Further dilutions should be made in a buffer containing a carrier protein.
<b>Preparation</b>	Recombinant protein expressed in <i>Pichia pastoris</i>
<b>Source</b>	<i>Pichia pastoris</i>
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Endotoxin Level</b>	< 1.1 EU/ug

**External Database  
Links**

**UniProt:**

[P51742](#)   [Related reagents](#)

**Entrez Gene:**

[403922](#)   TNF   [Related reagents](#)

---

**Synonyms**      TNFA, TNFSF2

---

**Product Information**      Canine tumor necrosis factor alpha (TNF alpha) is produced by a range of cells, including macrophages and monocytes and is an important mediator of the canine inflammatory response. Recombinant dog TNF alpha has an ED<sub>50</sub> of 10 pg/ml, determined by measuring the cytolytic effect of recombinant canine TNF alpha on murine WEHI-13VAR-164 cells in the presence of actinomycin D.

---

**Protein Molecular Weight**      17.3 kDa

---

**Purity**      >95% by SDS PAGE analysis

---

**References**      1. Grob, D. *et al.* (2021) Canine Angiostrongylus vasorum-Induced Early Innate Immune Reactions Based on NETs Formation and Canine Vascular Endothelial Cell Activation *In Vitro*. [Biology. 10 \(5\): 427.](#)

---

**Storage**      Prior to reconstitution store at +4°C.  
Following reconstitution store at -70°C.  
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**      12 months from date of despatch

---

**Health And Safety Information**      Material Safety Datasheet documentation #10302 available at:  
<https://www.bio-rad-antibodies.com/SDS/PDP001>  
10302

---

**Regulatory**      For research purposes only

---

**North & South America**      Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M433666:241118'

**Printed on 18 Nov 2024**