

Datasheet: PDP009

**BATCH NUMBER 171586**

<b>Description:</b>	RECOMBINANT DOG IFN ALPHA
<b>Name:</b>	IFN ALPHA
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	5 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Dog
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
<b>Source</b>	<i>Pichia pastoris</i>
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	10% Trehalose
<b>Product Information</b>	<b>Recombinant Dog interferon alpha (IFNα)</b> is a recombinant protein of approximately

19.0 kDa corresponding to a 164 amino acid sequence from native canine IFN $\alpha$ .

Canine IFN $\alpha$  is a type I interferon belonging to the alpha/beta interferon family. IFN $\alpha$  is produced by macrophages and involved in defence against viral infections.

---

<b>Protein Molecular Weight</b>	19.0 kDa
---------------------------------	----------

---

<b>Purity</b>	>95% by SDS PAGE analysis
---------------	---------------------------

---

<b>ELISA</b>	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with <a href="#">AHP2976</a> as a coating antibody and <a href="#">AHP2976B</a> as a detection antibody.
--------------	--

---

<b>Storage</b>	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.
----------------	--

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10591 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PDP009">https://www.bio-rad-antibodies.com/SDS/PDP009</a>
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M415625:230113'

Printed on 20 Oct 2025