

Datasheet: PDP014

BATCH NUMBER 170737

Description:	RECOMBINANT DOG INTERLEUKIN-1 BETA
Name:	IL-1 BETA
Other names:	IL-1F2
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Dog
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Source	<i>Pichia pastoris</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose

**External Database
Links**

UniProt:

[Q98292](#)

[Related reagents](#)

Product Information **Recombinant Dog interleukin-1 beta (IL-1 β)** is a recombinant protein corresponding to a 152 amino acid sequence of native canine IL-1 β .

Interleukin-1 β is a pleiotropic pro-inflammatory cytokine and member of the IL-1 superfamily. It is secreted by activated macrophages and monocytes and is involved in the modulation of a range of immune and inflammatory responses.

IL-1 β and the closely related IL-1 α , are both involved in haematopoiesis and exhibit many identical immunomodulatory properties, binding to the same type I interleukin-1 receptor, also known as CD121a.

Protein Molecular Weight 17.5 kDa

Purity >95% by SDS PAGE analysis

ELISA This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2980](#) as a coating antibody and [AHP2980B](#) as a detection antibody.

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10591 available at: <https://www.bio-rad-antibodies.com/SDS/PDP014>

Regulatory For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M415628:230113'

Printed on 20 Oct 2025