

Datasheet: PDP016

| | |
|----------------------|--------------------------------|
| Description: | RECOMBINANT DOG INTERLEUKIN-17 |
| Name: | IL-17A |
| Format: | Rec. Protein |
| Product Type: | Recombinant Protein |
| Quantity: | 5 µg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |
| Western Blotting | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|--|
| Target Species | Dog |
| Product Form | Purified recombinant protein - lyophilized |
| Reconstitution | Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. |
| Preparation | Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 10% Trehalose |
| External Database Links | UniProt: C6L8D7 Related reagents |

Product Information **Recombinant Dog interleukin-17A (IL-17A)** is a recombinant protein corresponding to a 127 amino acid sequence of native canine IL-17A.

Native canine IL-17A is a 17.7 kDa proinflammatory cytokine belonging to the IL-17 cytokine family and produced by activated T cells.

Protein Molecular Weight 14.8 kDa

Purity >95% by SDS PAGE analysis

ELISA This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2982](#) as a coating antibody and [AHP2982B](#) as a detection antibody.

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10591 available at:
10591: <https://www.bio-rad-antibodies.com/uploads/MSDS/10591.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M334910:181203'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)