

Datasheet: PDP020

BATCH NUMBER 153949

Description:	RECOMBINANT DOG INTERLEUKIN-13
Name:	IL-13
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Dog
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose

External Database Links

UniProt:

[Q9N0W9](http://www.uniprot.org/entry/Q9N0W9)

[Related reagents](#)

Entrez Gene:[442990](#) IL13 [Related reagents](#)

Product Information **Recombinant Dog Interleukin-13 (IL-13)** is a recombinant protein corresponding to aa 19-131 of native dog IL-13.

Dog interleukin-13 is an immunoregulatory cytokine produced by activated T lymphocytes. In dogs IL-13 is thought to have a regulatory role in allergic diseases. In in vitro experiments IL-13 was able to stimulate production of allergen-specific-IgE.

Protein Molecular Weight 12.3 kDa

Purity >95% by SDS PAGE analysis

ELISA This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2823](#) as a coating antibody and [AHP2823B](#) as a detection antibody.

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10591 available at:
<https://www.bio-rad-antibodies.com/SDS/PDP02010591>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373624:200923'

Printed on 18 Jan 2024