

Datasheet: PEP001

BATCH NUMBER 159837

| | |
|----------------------|--------------------|
| Description: | PURIFIED HORSE IgG |
| Name: | IgG |
| Format: | Purified |
| Product Type: | Purified Protein |
| Quantity: | 10 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-----------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |
| Immunodiffusion | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|--|
| Target Species | Horse |
| Product Form | Purified Equine IgG - liquid |
| Preparation | Purified IgG prepared by fractionation and DEAE chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃) |
| Approx. Protein Concentrations | IgG concentration 1.0 mg/ml |
| RRID | AB_324787 |

Product Information This product is a preparation of purified polyclonal **Horse IgG** prepared from normal horse serum.
By immunoelectrophoresis, this product was shown to react with antiserum specific for

horse IgGa and IgGb. No reactivity was observed for horse IgA, IgM or IgG(T).

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10258 available at:
<https://www.bio-rad-antibodies.com/SDS/PEP001>
10258

Regulatory For research purposes only

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M369477:200529'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)